

CSCI 14B

G3/14

BUNNYVALE, CA

FINAL REPORT
Section II
OPERATORS MANUAL
FOR
THE MAGNETOGRAPH PROGRAM

May 1974

Prepared for
George C. Marshall Space Flight Center
Huntsville, Alabama 35812

Contract No. NAS8-28018

Principal Investigator: Dr. A. M. Title

Lockheed Solar Observatory
Lockheed Palo Alto Research Laboratory
3251 Hanover Street
Palo Alto, California 94304

OPERATOR'S MANUAL FOR MAGNETOGRAPH PROGRAM

Larry November, August 21, 1973

To generate a directional magnetogram from original Spectra-Spectroheliograph data is a four pass process. Using the microdensitometer PASS1 scans the film to produce disk files numbered sequentially which contain the raw scan data. PASS2 may be used at any time in the process (i.e. when the disk gets full) to organize the data into a workable form and put this data on magtape. After all the data has been put on magtape PASS3 can be started which determines the field strength, angle, and velocities at each point on the sun. Near the end of the PASS3 run the computer requests a new magtape to output the final magnetograph. PASS4 can be used to examine the magnetograph final tape to produce any of eight types of magnetogram pictures (using the microdensitometer) or for types of vector magnetograms (using the Tektronix display scan).

This manual for use on the magnetograph program describes; 1. black box use of the program, 2. the magtape data formats used, 3. the adjustable control parameters in the program, and 4. the algorithms. With no adjustments on the control parameters this program may be used purely as a black box. For optimal use, however, the control parameters may be varied. The magtape data formats are of use in adopting other programs to look at raw data or final magnetograph data. For completeness I have included elaborate descriptions of how things work.

The sections of this manual are ordered as follows:

1. Using the Program - A Typical Run

- 1.1 About the Film - IX Coordinate and Polarization Pair
- 1.2 Densitometer Setup
- 1.3 Initializing the Computer System
- 1.4 A Warning
- 1.5 Zeroing the Microdensitometer
- 1.6 Labeling the Output File
- 1.7 User Defined Input Parameters
- 1.8 The Densitometer Scan
- 1.9 PASS2
- 1.10 The Next Densitometer Scan
- 1.11 After the Data Tape is Done
- 1.12 PASS4

2. Data Formats

- 2.1 Disk Data Files after PASS1
- 2.2 Data Tape Format after PASS2
- 2.3 Data Tape Format for Magnetogram after PASS 3
- 2.4 Disk Picture Files from PASS4

3. Control Parameters and Program Relinking

- 3.1 PASS2C.FTN
- 3.2 PASS3D.FTN

4. Algorithms

4.1 PASS2

4.2 PASS3

4.3 PASS4

1. Using the Program - A Typical Run

Given spectra-spectroheliograph data covering a region on the sun, information is available to readily determine velocities, B field strengths, and B field angles. This section describes a four pass computer process which is capable of deriving this information from the raw data. The sub-sections of this part of the manual are sequentially ordered to present the complete process of how to make a final vector magnetogram from the raw film data.

1.1 About the Film - IX Coordinate and Polarization Pair

The spectra-spectroheliograph film has two basic formats, the old type, and the new type, which are shown in Figure 1.1. The basic characteristics of these two types are different and it is important when running PASS1 to distinguish between them.

A COMPLETE RUN is composed of several REGION SCANS which are each composed of many SLIT PAIRS. Before each REGION SCAN the spectra-spectroheliograph was pointed at some initial solar coordinate and allowed to step across the region making many SLIT PAIRS each at a different X coordinate. After a change in polaroids the next REGION SCAN was begun. In each REGION SCAN the same pattern of spots was crossed and a similar looking sequence of SLIT PAIRS encountered.

It is extremely important that the film be correctly identified or the analysis cannot correctly proceed. Each SLIT PAIR is identified by an X coordinate integer which I refer to as IX. IX sequentially increases as each SLIT PAIR in a REGION SCAN is crossed. Also IX is the same for SLIT PAIRS that were taken at the same solar coordinate in separate REGION SCAN's of the same COMPLETE RUN. It is essential that (1) IX increase sequentially and (2) that an accurate correspondence exist between the SLIT PAIR numbering and the actual solar coordinates.

Sometimes only a portion of a REGION SCAN needs be analyzed which places one more restriction on IX. I introduce the term ANALYZING SCAN to represent that portion of the REGION SCAN that is analyzed. Just as there may be four REGION SCANS for a COMPLETE RUN there may be (at most) four ANALYZING SCANS. Figure 1.2 portrays how the ANALYZING SCAN fits in. Note that each ANALYZING SCAN in a COMPLETE RUN contains the same set of IX coordinates. Also note that the first SCAN PAIR of the ANALYZING SCAN is labeled IX=1. This is the third and last restriction on IX.

In addition to IX each SCAN PAIR has another characteristic associated with it. This characteristic is its polarization. Each pair physically represents the light in two mutually orthogonal polarizations. Thus a SCAN PAIR may represent RHC-LHC polarizations, or linear 37° - linear 127° , a linear 0° - linear -90° . Each REGION SCAN has associated with it one such polarization pair, and all the SCAN PAIRS in that REGION SCAN WERE TAKEN through the same pair of polaroids.

Each SCAN PAIR can thus be uniquely identified by IX and its polarization pair.

1.2 Densitometer Setup

Densitometer Setup requires two operations. First the film must be properly cleaned and placed on the platten. Second the microdensitometer must be calibrated and focused for the film.

The film should be cleaned and placed on the platten such that the REGION SCAN begins on the right and ends somewhere in the roll on the left. IX = 1 should be placed somewhere near the right edge of the platten with IX increasing to the left (as the densitometer is faced). The film must be aligned such that the film edge lies precisely parallel to the microdensitometer X axis. This can be checked by running the densitometer back and forth in manual.

Once aligned the densitometer may be precisely focused, a small aperture installed and aligned, and the transmission set for a clear portion of the film. The exact procedure is not described here as an elaborate method is given in the 'Microdensitometer System' by Steve Schoolman.

1.3 Initializing the Computer System

Using the program PIP on the PDP-11 system initialization of the DISK and the MAGTAPE can be accomplished. A SCRATCH disk should be installed and running on DKØ, and a blank tape running on MTØ.

Five restrictions exist for the scratch disk. First, at least 3000 blocks of free disk must be freed. Second, a file called NUMBER must be deleted from the user's directory. Third, a file called PARM must be deleted from the user's directory. All these conditions can be met using PIP. Third, the file PASS1.LDA must be present on the disk. Fourth, the file PASS2.LDA must be present on the disk. Both of these files can be found on MAGTAPE S-24 under the User Initials [2,11] (please do not zero this tape).

It is necessary that a ZERO mark be present on the output tape. If one is, no more action need be taken. If one is not the command MT0=/ZE in PIP will zero the tape (and wipe out anything already on the tape).

Once both the disk and tape are readied the program is ready to begin.

1.4 A Warning

The program PASS2 has defined within a BLOCK DATA statement all the program control parameters. In general those parameters should work for any piece of film. However, a new dimension in film size or raw step wedge data could lead to trouble. It may be necessary under unusual conditions to alter the control parameters and re-link PASS2 Section 3 of this booklet describes the conditions and the procedures under which PASS2 should be altered.

1.5 Zeroing the Microdensitometer

Zeroing is accomplished with the program PASS1.LDA which should be called in from disk and running. PASS1 is just the system program TRACE with one alteration - it smooths data as it scans which TRACE does not do. PASS1 (and TRACE) is a user interaction program with many features only a few of which I describe here. PASS1 is in idle and waiting for a user response if a * is the last display character.

With the densitometer in MANUEL the densitometer γ should be moved to STARTING POSITION 1 on the film (see Figure 1.2). Z is the command which will then zero the reading to this starting point.

1.6 Labeling the Output File

The command I in PASS1 causes the computer to respond with:

IDENT:

After which the user can enter a label of 31 characters length or less. A typical label identifies the region and the time of the COMPLETE RUN.

1.7 User Defined Input Parameters

Essentially three decisions are now left up to the user. He may select the parameter L = length of film he wants to scan (microns), ΔX = the resolution

needed in X (microns), and ΔY = the resolution reading in Y (microns). Typically these parameters are on the order of $L = 15\text{cm} = 150,000 \mu$, $\Delta X = 15\mu$, $\Delta Y = 200\mu$. Having decided on these parameters the user-computer dialogue should proceed as follows (computer response is underlined):

For old data:

* U

USER-DEFINED SCAN PARAMETERS:

X-DIR= +

Y-DIR= -

PATTERN= R or E (user's choice)

DELTA X= value of ΔX <CR>

POINTS/LINE= value of $L/\Delta X$ ($\leq 24,576$) <CR>

YSTEP= value of ΔY <CR>

LINES= value of $57150/\Delta Y$ <CR>

X= 0 <CR>

Y= 0 <CR>

(remember the line feed termination character)

For new data:

* U

USER-DEFINED SCAN PARAMETERS:

X-DIR= +

Y-DIR= -

PATTERN= R or E (user's choice)

DELTA X= value of ΔX <CR>
POINTS/LINE= value of $L/\Delta X$ (24,576) <CR>
YSTEP= value of ΔY <CR>
LINES= value of $30,000./\Delta Y$ <CR>
X= \emptyset <CR>
Y= \emptyset <CR>
X= \emptyset <CR>
Y= -30480 <LF>
 (remember the line feed).

If errors are made on this input it may be corrected by typing U again - this time the computer will respond:

* U

OPTION:

Correction may be made to any of the parameters by selecting the appropriate option. These are the options:

U - complete dialogue
 X - X direction
 Y - Y direction
 P - Pattern
 C - Starting coordinates
 DX - Delta X
 DY - Y step
 NP - Number of points per line
 NL - Number of lines
 E - exit to monitor from PASS1.

1.8 The DENSITOMETER SCAN

Now the actual DENSITOMETER SCAN can begin. The initiating character for the scan is S. PASS1 first opens a file for output of the data on disk, and then proceeds to scan recording the data on disk. At the end of the scan a * flashes on the display screen and the process may continue.

1.9 PASS2

Generally speaking one densitometer scan is enough to fill up the disk. PASS2 should therefore be called after each densitometer scan to move the data onto tape.

PASS2 first will ask to know which files it is to work on off of disk. This is specified by the letter which identifies the file on disk and the number of the first and last file to be analyzed from disk (these two numbers may be equal if only one file is to be analyzed).

When PASS2 is brought into core for the first time during a COMPLETE RUN it must be initialized with certain parameters. It will ask for the following parameters:

WAVELENGTH SCALE (MICRONS WAVELENGTH/MICRON ON FILM)

($\approx 1.40000E-6$)

XSTEP (SECONDS OF ARC /SLIT PAIR) (≈ 0.5)

YUNIT (SECONDS OF ARC/MICRON ON FILM) ($\approx 2.4038 \times 10^{-3}$)
WAVELENGTH OF FILM LEADING EDGE (MICRONS) (≈ 52.4994)
WAVELENGTH OF LINE CENTER (MICRONS) (≈ 52.502)
WAVELENGTH OF ANY CALIBRATION LINES (NONE - <CR>)

Next PASS2 will tell which data file it has opened and request input of the initial IX value for the densitometer scan. IX=1 for the first scan or in general it is equal to the value of IX for the first SLIT PAIR in the DENSITOMETER SCAN.

PASS2 will next request the two orthogonal polarizations of the SLIT PAIRS in this DENSITOMETER SCAN. The information for RHC polarization is coded 511, and for LHC polarization 256. Linear polarization is entered by the angle of the polaroid in degrees from the X-scan direction on the sun.

PASS2 will then proceed to move this data onto tape in a more convenient format. As a warning PASS2 requires 100 contiguous scratch disk blocks. If it cannot allocate this space a STOP ~~000001~~ will occur. This situation can be remedied by using the /PK switch in PIP or by deleting more files on disk and rerunning PASS2. PASS2 takes approximately 15 minutes to run through a full DENSITOMETER SCAN file.

1.10 The Next DENSITOMETER SCAN

Once PASS2 and PASS1 have run the files PARM and NUMBER should be present on disk. Do not delete either of these two files until the end of the COMPLETE RUN. Once again, though, room must be made on disk which can be accomplished by deleting the old data file. Once the disk is readied the next DENSITOMETER SCAN may begin.

Move the film on the PLATTEN to the next starting point on the film (see Figure 1.2). Once again the Densitometer must be properly calibrated and the film set parallel to the Densitometer X axis (see Section 1.2). Once done, PASS1 must be rerun and zeroing (section 1.5), labeling (section 1.6), and initialization of input parameters (section 1.7) redone. Please note that the input parameters L and ΔX may vary from DENSITOMETER SCAN to DENSITOMETER SCAN, but that ΔY must not change. (as L, ΔX , ΔY are defined in section 1.7). After all is readied S (see section 1.8) initializes the DENSITOMETER SCAN.

After the SCAN PASS2 should be run. This time through, PASS2 will ask only for the data file letter and number, IX for the first SLIT PAIR in this DENSITOMETER SCAN, and for the orthogonal polarizations of this DENSITOMETER SCAN. After PASS2 finishes this DENSITOMETER SCAN will have been completely moved to tape and the next DENSITOMETER SCAN begun.

1.11 After the Data Tape is Done

After the last of the DENSITOMETER SCANS in the COMPLETE RUN has been moved to tape by PASS2 the actual magnetic field information can be processed. PASS3 performs this function.

PASS3 is recorded on tape S-24 under the UIC (2,11) and must be moved to disk before being run. Three considerations should be made before running PASS3. First, for optimal running the control parameters should be correctly adjusted (see section 3). Second, PASS3 requires (number of points in X) * number of points in Y) / 64 blocks of contiguous scratch area on disk which must be available. Third a tape with a correct zero mark should be ready and waiting for output of the final magnetogram.

When the disk is ready and the data tape is mounted on MTØ and moved to the load point then PASS3 may be run. PASS3 first requests the maximum picture dimensions and then a single letter which defines the files on the data tape to be processed. PASS3 then proceeds to process the data which takes approximately (Number points in X) * (Number of Points in Y) / 200 minutes.

At the completion of the data processing, PASS3 will ring a bell and request the final magnetogram tape be mounted on unit MTØ. Once loaded and an output file is assigned to unit 8, PASS3 should be continued. After output of the magnetogram PASS3 will rewind the tape and return control to the monitor.

1.12 PASS4

Final display of the data from the magnetogram tape in various formats is possible using PASS4. PASS4 requires only that the magnetogram tape be mounted on MT0 and that sufficient space is available on disk for output of the picture.

When PASS4 is run it first will request assignment of the magnetogram file to tape unit 8. After a continuation is given by the user PASS4 will request two integers, MODE and ITYPE. MODE=0 causes an exit to monitor. The other MODE and ITYPE possibilities are described:

MODE=1, ITYPE=1 - A dopplergram will be produced. The dopplergram is written on a disk file whose name is given by the program. This file may later be output on the densitometer as a picture.

MODE=1, ITYPE=2 - A magnetogram will be produced. The magnetogram is written on a disk file which may later be rendered as a photograph by the densitometer.

MODE=1, ITYPE=3 - A gamma-gram is produced. This disk file may later be rendered on the densitometer as a photograph which displays fields toward the observer as black (0° yields transmission of 0) and fields away from the observer as white (180° yields transmission of 1023).

MODE=1, ITYPE=4 - A phi-gram is produced. This disk file may later be rendered on the densitometer as a photograph which displays fields in the X scan direction as black (0° yields transmission of 0) and fields in the -X scan direction as white (180° yields transmission of 180°).

MODE =2, ITYPE=1 - A longitudinal magnetogram is produced. $B \cos (\delta)$ is written on the output file which may later be rendered as a photograph.

MODE =2, ITYPE=2 - A transverse magnetogram is produced. $B \sin (\delta)$ is written on the output file which may later be rendered as a photograph.

MODE =3, ITYPE=1 - Three longitudinal magnetograms are produced each representing a different color. $B \cos (\delta)$ is written on the output files, and the color wheel red-blue-yellow represents the field angle of going from 0° to 180° from the X scan direction.

MODE=3, ITYPE=2 - Three transverse magnetograms are produced each representing a different color. $B \sin (\delta)$ is written on the output files, and the color wheel red-blue-yellow represents the field angle δ going from 0° to 180° from the X scan direction.

MODE=4, ITYPE=1 - A vector magnetogram is displayed on the tektronix screen as seen from the viewing angle (δ, δ') as requested by the program. All fields are displayed which are greater than B_{MIN} (program request) as vectors of equal length which point in the direction of the field and are anchored to the surface of the sun in the appropriate location. After completion of the drawing the bell is sounded. Striking a return on the keyboard will cause a hardcopy to be produced and the program to continue.

MODE=5, ITYPE=1 - A vector magnetogram is displayed on the tektronix screen as seen from the viewing angle (δ, δ') (program request). All fields are displayed which are greater than B_{MIN} (program request) as vectors of relative length

which point in the direction of the field and are anchored to the surface of the sun in the appropriate location. After completion of the drawing the bell is sounded and a return from the keyboard will cause a hardcopy to be produced.

The files that have been created by PASS⁴ in modes 1 through 3 may be rendered as photographs by the densitometer using the system program TRACE. Setup of the densitometer is described in the "Operator's Manual for Microdensitometer Control Program" by Steve Schoolman. The film should be on the platten with the densitometer at the lower-left of the film. Playback of the photograph is initiated by the P command in running TRACE.

\$ RUN TRACE

MONITOR

* Z

* P

the computer then responds and the dialogue continues:

SCALE FACTOR user's choice 1 to 100 <CR>

PARAMETER SOURCE? R

FILE NAME: complete name of file created by PASS⁴ <CR>

TYPE ANY KEY TO CONTINUE

With the lights out and the film placed on the platten the user may initiate the picture output by typing any character. Completion of the picture output is signified by a bell.

2. Data Formats

Often times the user may want to use the output data or the output magnetogram in ways that are not possible through this system. In order to allow for easy use of the output information I now describe the data formats which are used.

I describe here a total of four data formats. They are: 1) The format of the data at the end of PASS1 as it is stored on disk, 2) the format of the data tape as it is produced by PASS2, 3) the format of the magnetogram tape as produced by PASS3, and 4) the format of the picture records as produced by PASS4 in modes 1 through 3.

Each of these total mediums will contain one or many files each of which will contain many Fortran unformatted records. Each record must be read in a single Fortran READ statement and may contain as many as 12000 integer words (2 bytes each). All the records, in all the files, in all four storage media are of the same general format. The first word of each record is a code word, which describes the nature of the record, and the second word is a length word which tells how many words remain in the record. The rest of the record contains the desired information.

2.1 Disk Data Files after PASS1

Each disk data file contains microdensitometer data from one DENSITOMETER SCANS. A COMPLETE RUN is liable to consist of many DENSITOMETER SCANS and so many of these data files.

The data file consists of many records. The first two or three of these records contains control information and then each additional record contains the raw densitometer data from each line of the scan, one record to each line. The record contains the following.

RECORD 1 = Label:

WORD 1 = ICODE = 1 or 4

WORD 2 = N

WORD 3 to WORD N+2 = a label 2N bytes long

RECORD 2 = Message (optional)

WORD 1 = ICODE = 2

WORD 2 = N

WORD 3 to WORD N+2 = a message 2N bytes long

RECORD 2 or 3 = Parameters

WORD 1 = ICODE = 3

WORD 2 = N = 141

WORD 3 = NPPL = Number of points per DENSITOMETER SCAN line.

WORD 6 = IDELX = distance between points in microns.

WORD 7 = IDELY = distance between lines in microns.

WORD 10 = NLPF = number of lines per frame (for old data this is the number of lines for new data this is half the number of lines).

WORD 11 = -1 for raster pattern, 0 for edge scan pattern.

WORD 14 = 0 for old data, 2 for new data.

RECORD 3+ or 4+ = Densitometer Data.

WORD 1 = ICODE = - line number

WORD 2 = N = NPPL

WORD 3 to WORD N+2 = the data.

More information is available in the parameter record than is defined here.

A complete description can be found in the output Format section of the 'Operator's Manual to Trace' by S. Schoolman.

2.2 Data Tape Format After PASS2

The data tape produced by PASS2 contains a sequence of files labeled <letter> n.DAT where n varies from 1 to the total number of DENSITOMETER SCANS in the COMPLETE RUN. Each of these files contains the following sequence of records:

RECORD 1 = Label

WORD 1 = ICODE =1

WORD 2 = N

WORD 3 to WORD N+2 = label of length 2N bytes.

RECORD 2 = Parameters

WORD 1 = ICODE=2

WORD 2 = N = 10

WORD 3 and 4 = XSTEP = seconds of arc per step in scan direction
real number four bytes long).

RECORD 2 = Parameters

WORD 1 = ICODE = 2

WORD 2 = N = 6

WORD 3 & 4 = XSTEP = seconds of arc per step in X (real, four bytes long).

WORD 5 & 6 = YSTEP = seconds of arc per step in Y (real).

WORD 7 = NX = total number of steps in X direction.

WORD 8 = NY = total number of steps in Y direction.

RECORD 3+ = Data

WORD 1 = ICODE = 0.

WORD 2 = N = 10

WORD 3 = IX = number of X step for this data point.

WORD 4 = IY = number of Y step for this data point.

WORD 5 & 6 = V = velocity (Km/sec, positive-downward) (real).

WORD 7 & 8 = B = field strength (gauss) (real).

WORD 9 & 10 = GAM = field angle in degrees from the line of sight (real).

WORD 11 & 12 = PHI = field angle in degrees from X (real).

2.4 Disk Picture Files from PASS4

Each picture is encoded on one file. The records are in the same format as those output by PASS1, i.e. capable of controlling the microdensitometer program

TRACE:

RECORD 1 = Label

WORD 1 = ICODE = 4

WORD 2 = N

WORD 3 to WORD N+2 = label of length 2N bytes.

RECORD 2 = Message

WORD 1 = ICODE = 2

WORD 2 = N

WORD 3 to WORD N+2 = message.

RECORD 3 = Parameters

WORD 1 = ICODE = 3

WORD 2 = N = 141

WORD 3 = NPPL = number of points per line.

WORD 6 = IDELX = number of microns between output points in X=10.

WORD 7 = IDELY = number of microns between output points in Y=10.

WORD 10 = NL = number of lines in output picture.

WORD 12 = ISPEED = 64

All other words in this record are zero.

RECORD 4+ = Data

WORD 1 = ICODE = -line number

WORD 2 = N

WORD 3 to WORD N+2 = transmissions of points on line.

3. Control Parameters and Program Relinking.

On tape S-24 all the necessary programs are stored for making a magnetogram.

S-24 contains the following files under UIC [2,11]:

PASS1.LDA	PASS3A.FTN	DELETE.MAC
PASS2.LDA	PASS3B.FTN	
PASS3.LDA	PASS3C.FTN	
PASS4.LDA	PASS3D.FTN	
PASS2A.FTN	PASS4A.FTN	
PASS2B.FTN	PASS4B.FTN	
PASS2C.FTN	PASS4C.FTN	

Direct running can be accomplished using the default control parameters in the load modules. The control parameters for PASS2 are contained in PASS2C.FTN and the control parameters for PASS3 are contained in PASS3D.FTN. If these files are altered then PASS2 or PASS3 must be recompiled with the /ON switch and relinked. If PASS2C.FTN is changed PASS2 must be linked from PASSA/CC, PASS2B/CC, PASS 2E, and DELETE. If Pass 3D.FTN is altered then PASS 3 must be relinked from PASS3A/CC, PASS3B/CC, PASS3C/CC, PASS3D, DELETE.

Listings of these programs can be found at the end of this booklet. The listings are by file and give all the subprograms and programs used.

3.1 PASS2C.FTN

PASS2C is a BLOCK DATA statement which defines parameters in four common blocks.

Each of these four commons is described with the parameters:

/FILM/ contains all the information on the film dimensions. In practice this block is the most critical and these parameters must be changed in switching program use from old to new data. Five parameters in /FILM/ are adjustable SLT(1), SLT(2), SMARG(1), SMARG(2), SEFF. Figure 3 illustrates these dimensions for the new and the old film. Each parameter is described below:

SLT(1)= size in microns of film of one slit crossing (default= old data = 2644.0)

SLT(2) = size in microns of film of second slit crossing in SLIT PAIR (default = old data = 2644.0)

SMARG(1) = distance in microns of film between second slit crossing in SLIT PAIR to first slit crossing of next pair (default = old data = 3527.0)

SMARG(2) = distance in microns of film between first slit crossing in SLIT PAIR to second slit crossing in next pair (default = old data = 685)

SEFF = distance around line center in microns wavelength which is useful - anything over this distance will be discoded (default = 1)

As a check 9500.0 microns represents the total recycling distance. No matter what value for SLT(1 & 2) and SMARG (1 & 2) are used this check applies.

On the new film format these four parameters have the values:

SLT(1) = 8500.

SLT(2) = 8500.

SMARG(1) = 1000.

SMARG(2) = 1000.

/DENS/ contains only one parameter ICONTR. The program in attempting to find the leading edge of the film looks for a slope in transmission of ICONTR over ten points. For old data it looks for a drop in transmission of ICONTR; for new data, it looks for a rise in transmission of this amount. If a STOP2 of the program results it is likely that this parameter was too large.

/WEDGE/ contains step wedge data if it is available. Two integer arrays are defined in /WEDGE/, ITRANS(11) and INT(11). ITRANS(n) represents the transmission reading that results for the corrected transmission INT(n). These two relations should be approximately linear, i.e. INT is approximately proportional to TRANS and ranges from about 0 to about 1000. NSTEP tells the number of useful steps in the wedge. If no correction is necessary or no wedge data is available then NOWEG = .TRUE. The default for this common is NOWEG = .TRUE.

/BUFFER/ has only one adjustable number - that is the size of the array IBUF. This number represents the maximum densitometer line length and is determined by the amount of core available on the machine. The default value of 12288 is about the maximum that will run in a 24K core machine.

3.2 PASS3D.FIN

PASS3D is a BLOCK DATA statement that contains program control parameters in four commons. These four commons are /PROFILE/, /LNPRM/, /FILTER/, /FFT/ and are described as follows:

/PROFILE/ contains two numbers which characterize the line profile. XWIDTH represents the half-height field width of the line in microns wavelength.

PKLIN represents the relative line depth ($\sim 25 \times (\text{continuum intensity})/(\text{line center intensity})$). The intensity profile is assumed to be approximately lorentzian with the following relative curve:

$$I = \frac{PKLIN}{4 \left(\frac{\Delta \lambda}{XWIDTH} \right)^2 + 1}$$

The default values which the program uses are XWIDTH=.06E-2 microns wavelength and PKLIN=48.

/LNPRM/ contains other line parameters that apply to splitting and the doppler line shift. Only one control parameter is settable in this block; that is GFAC = Lande' g factor for the line in question. A default value of GFAC=3. is assumed by the program.

The program is designed to act like many filters which sample the data at many points. /FILTER/ controls the nature of the filters. NSAMP, the first element

of /FILTER/ equals the number of bandpasses to be used by the program (≤ 15 and odd). XBAND is the bandpass of each filter in microns wavelength and XSTEP is the distance between bandpasses in microns wavelength. Default values are NSAMP=15, XBAND=.01E-2, XSTEP=.03E-2.

It is possible in this program to invoke an alternate routine for determining field strength. /FFT/ contains the control data for this routine. NOFFT=.TRUE. causes this technique to be skipped altogether. BMIN sets the lower limit of fields that this routine will consider. GBND gives the number of degrees around $\delta = 90^\circ$ in which this routine is not called. The default values are presently set at NOFFT=.TRUE., BMIN=1800. (gauss), GBND = 20.(degrees).

4. Algorithms

The details of operation of PASS2, PASS3, and PASS4 are described here. PASS2 and PASS4 are briefly and basically described. PASS 3 is described in greater rigor since this is where the real field determination goes on.

4.1 PASS2

PASS2 has the function of converting microdensitometer raw data into manageable data records. The original microdensitometer data consists of records that represent DENSITOMETER SCAN lines. The final output tape saves only the useful data in records, one slit crossing to a record.

PASS2 has the goal to keep track of label records, correctly set and use control parameters, and breakdown of the data. In PASS2A.FTN the MAIN program is listed. Initialization and the keeping track of control parameters is done by the first part of the program (up to statement 40). The loop DO 60 contains the data examination part of the program.

In the initialization phase two possible directions are taken. If this is the first time PASS2 is being run in a complete run (IFIL1=1) then the program requests initial parameters which it stores in a data file SPARM on disk. If this is not a first run (IFIL1/1) then PASS2 sets the initial parameters by checking the file SPARM already on disk.

The data reduction phase of loop DO 60 first opens data files for reading, then reads sequentially all the records on each of these files. As each record it reads it is processed according to what its code word ICODE implies. The actual breakdown of line data records is called if ICODE is negative. Under this condition the loop DO 54 reads the data from the raw data file in a random access technique (using the subroutine READ). PROCSS then handles each line as it is given to it.

PROCSS goes through a four step operation to break down a data line. First if the line was scanned backwards in a raster scan - the line is flipped (subroutine FLIP). Second the edge of the first slit crossing on the line is found by CR3DEG. Third, the line is stepped through and broken down according to the distance given in the common /FILM/ and results are written on the scratch file

SCRATCH.TMP. Fourth, after information for the two orthogonal polarizations is present in SCRATCH.TMP (takes two line process for new data and a one line process for old data) it is rewritten on tape in a finalized form.

This process continues until all the raw data files on disk are read. PASS2 then closes all its files and deletes SCRATCH.TMP.

4.2 PASS3

The heart of the operation is here. PASS3 takes the data tape created by PASS2 and yields magnetic field strengths, angles, and fluid velocities. In this section the mechanics of the operation is discussed. The theory behind the operation is described in the Fourier Transform method by A. Title and T. Tarbell.

Because of the severe space limitation on disk PASS3 must operate as a two step program. First data is read off of tape. As each data record is read all the essentials of that record (field, angles, velocity) are derived in the routine PROCSS using the routine PRMFFT. As results come out they are recorded in an extremely condensed matter on a large contiguous scratch file on disk. Finally after all the raw data has been looked at the computer requests a new tape be mounted for output of the final magnetogram. PASS3 then begins its second phase of operation and moves data from the SCRATCH.TMP file onto tape. Before exit SCRATCH.TMP is deleted.

It is vital to understand the structure of the in between scratch file to understand the program mechanics. Because it is contiguous this scratch file can be treated as a large array. IO1 is a routine which gives SCRATCH.TMP this appearance. IO1 references V, B, δ , ϕ according to the coordinate IX, IY. Thus any velocity may be easily read or written for any point in the picture at any time in the execution. As field values are found from the data they are averaged with previous results (if there are any and weighted properly), and returned to SCRATCH.TMP over the old location. During a COMPLETE RUN as many as four values for the field may be determined for each point in the picture. These have been appropriately handled so as to take only one word on disk during the COMPLETE RUN.

The work of PASS3 is carried out in PRMFFT and the routines it calls CONTRA, COMBINE, POLYAN, and BGOOD. In addition there are a set of operational routines that perform mathematical operation for PRMFFT, these are DSUM, COSTRA, ZER, PLINT.

Basically PRMFFT manages the fourier transform method of obtaining $|B|$, cosine γ , and the line of sight velocity. Subroutine CONTRA determines the continuum values. COMBINE determines the line center position, and hence the line of sight velocity; while POLYAN manages the determination of the zero crossing, and cosine γ . The value of the zero crossing is determined by ZER and COSTR calculates the required fourier transforms. The polynomial integration procedure for cosine gamma is performed in PLINT. Finally the values determined by the above routines are reviewed in BGOOD. The routine DSUM adds sets of values for a number of averaging processes.

4.3 PASS4

There are no tricks in PASS4. The algorithms are obvious and straightforward. Given magnetogram data PASS4 can be used to construct magnetogram pictures in any of twelve formats.

MAIN (PASS4A) serves the function of determining which output feature is desired by the user. MAIN only sets certain variables which are then passed to the appropriate routines for handling.

PICTUR creates a file for the densitometer which can later be interpreted as a photograph. It uses the two routines LABEL, and SETSCL to write the label and parameter block of the output file. LNINT is where the work goes on, as this routine determines the actual intensity for the picture point being written. PICTUR may be called when MODE 1, 2, or 3 of the routine are entered.

MODES 4 and 5 of the program cause the routine VECT. to be called from MAIN. VECT examines the magnetogram data and constructs on the output screen a vector magnetogram according to specifications.





Figure 1.1 Old and New Film Formats for Spectra-Spectraheliograms.



This is an example of a piece of spectra-spectraheliograph film with the terms identified. In the complete run of this example the region was scanned three times in three polarization pairs. REGION SCAN 1 was LHC, RHC. (IPOL = 256, 511), REGION SCAN 2 was x, y (IPOL = 0, 90), REGION SCAN 3 was in linear polarizations 45° , 135° (IPOL = 45, 135). Corresponding to the three REGION SCANS are three ANALYZING SCANS each of which begin at the same point on the sun. This solar starting point is labeled IX = 1. In this example each ANALYZING SCAN took three separate DENSITOMETER SCANS. Each DENSITOMETER SCAN covered about 7 slit pairs (not clear in drawing, but may vary even in real run).




Figure 1.2

```

C PASS 2 -DATA FROM DISK TO TAPE
C THIS PROGRAM TAKES A SEQUENCE OF M INPUT FILES LABELED <CHAR>1.DAT
C TO <CHAR>M.DAT OF DENSITIZED DATA AND PRODUCES AN OUTPUT TAPE
C OF IDENTICALLY LABELED FILES WITH THE DATA ORGANIZED INTO SMALL
C RECORDS EACH CONTAINING ONE SLIT CROSSING. THE OUTPUT DATA TAPE
C CONTAINS SEQUENTIAL FILES EACH FILE CONTAINING MANY RECORDS.
C EACH RECORD IS IDENTIFIED BY ITS FIRST WORD =ICODE AND BY ITS
C SECOND WORD N WHICH TELLS THE LENGTH OF THE BUFFER TO
C FOLLOW. WITHIN THE BUFFER THE INFORMATION IS CODED AS FOLLOWS:
C ICODE=1,N,IBUF=LABEL OF LENGTH N WORDS.
C ICODE=2,N,IBUF(1 AND 2)=XSTEP=SECONDS OF ARC PER X STEP(REAL).
C IBUF(3 AND 4)=YSTEP=SECONDS OF ARC PER Y STEP.
C IBUF(5 AND 6)=ZSTEP=MICRONS OF WAVELENGTH PER DATA POINT.
C IBUF(7 AND 8)=ZCEN=MICRONS WAVELENGTH OF CENTERLINE.
C IBUF(9)=NLPF=NUMBER OF YSTEPS IN Y DIMENSION.
C IBUF(10)=1 IF OLD DATA, 2 IF NEW DATA.
C ICODE=0,N,IBUF(1)=J=NUMBER OF POINTS IN DATA RECORD=N*6.
C IBUF(2)=IX=NUMBER OF X STEP.
C IBUF(3)=IY=NUMBER OF Y STEP.
C IBUF(4 AND 5)=Z0=NUMBER OF Z STEPS (REAL) FROM FIRST RECORD TO LINE
C CENTER.
C IBUF(6)=IPOL=POLARIZATION OF RECORD (511=RHC,256=LHC,0 TO 180=
C LINEAR).
C IBUF(7) TO IBUF(N)=FILM DENSITIES, SMOOTHED AND PHOTOMETERED.
C RECORDS ARE ORDERED IN ORTHOGONAL POLARITY PAIRS SEQUENTIALLY
C IN IX FOR K STEPS, FOR EACH VALUE OF IY.
C LARRY NOVEMBER = 21-AUG-73
COMMON /BUFFER/ICODE,IBUF(1)
DIMENSION ZLNE(5),IPOL(2),NUM(8),IEXT(5)
LOGICAL END
LOGICAL*1 IFIL(10),IFILNM(4),CHAR(2)
EQUIVALENCE (IFIL(1),IEXT(1)),(IFIL(3),IFILNM(1)),
* (CHAR(1),ICHAR)
DATA IEXT/'00','00','00','0.D','AT'/,ICHAR/' 0'/
CALL SETFIL(2,'SCRATCH.TMP',IERR,'DK',0,0,'233,2)
DEFINE FILE 2(50,512,U,IVAR)
WRITE(2,1,ERR=3) IDUMMY
GO TO 4
3 CONTINUE
END FILE 2
WRITE(6,1000) IERR
STOP 1
4 CONTINUE
WRITE(6,1000)
READ(6,902) IFIL(1),IFIL1,IFIL2
IF(IFIL1.EQ.1) GO TO 6
CALL SETFIL(1,'SPARM')
READ(1) NLNE,ZUNIT,XSTEP,YUNIT,ZLNE
END FILE 1
GO TO 12
6 CONTINUE
WRITE(6,1001)
READ(6,900) ZUNIT,XSTEP,YUNIT
WRITE(6,1002)
READ(6,900) ZLNE(1),ZLNE(2)

```

```

0028      NLNE=0
0029      DO 10 L1=3,5
0030      WRITE(6,1003)
0031      READ(6,900) ZLNE(L1)
0032      IF(ZLNE(L1).EQ.0.) GO TO 11
0033      NLNE=NLNE+1
0034  10 CONTINUE
0035  11 CALL SETFIL(1,'SPARM')
0036      WRITE(1) NLNE,ZUNIT,XSTEP,YUNIT,ZLNE
0037      END FILE 1
0038  12 CONTINUE
0039      DO 60 L0=IFIL1,IFIL2
0040      ENCODE(4,1006,IFILNM) L0
0041      DO 13 L1=1,4
0042  13 IF(IFILNM(L1).EQ.CHAR(1)) IFILNM(L1)=CHAR(2)
0043      WRITE(6,1007) IFIL
0044      CALL SETFIL(1,IFIL)
0045      CALL SETFIL(8,IFIL)
0046      WRITE(6,1005)
0047      READ(6,901) IX
0048      WRITE(6,1004)
0049      READ(6,901) IPOL
0050      CALL CURVE
0051      IREC1=0
0052      IFR=0
0053  15 CONTINUE
0054      IREC1=IREC1+1
0055      CALL READ(IREC1,N,ICODE,IBUF,END)
0056      IF(END) GO TO 55
0057      IF(ICODE.LE.0) GO TO 51
0058      GO TO (30,50,50,30),ICODE
0059  30 ICODE=1
0060      WRITE(8) ICODE,N,(IBUF(L1),L1=1,N)
0061      GO TO 15
0062  50 CALL SPARM(IBUF,ZUNIT,XSTEP,YUNIT,NPPL,NLPF,NFR,
*      LPAT,ZLNE(2),ZSTEP)
0063      GO TO 15
0064  51 DO 54 L1=1,NLPF
0065      DO 54 L2=1,NFR
0066      CALL PROCSS(ICODE,IBUF,IPOL,IX,L1,ZSTEP,NLNE,ZLNE,
*      LPAT,NPPL,L2,NFR)
0067      CALL READ((NFR-L2)*NLPF+L1+IREC1,N,ICODE,IBUF,END)
0068      IF(END) GO TO 55
0069  54 CONTINUE
0070  55 CONTINUE
0071      END FILE 1
0072      END FILE 8
0073  60 CONTINUE
0074      END FILE 2
0075      CALL DELETE
0076  900 FORMAT(E10.0)
0077  901 FORMAT(I3)
0078  902 FORMAT(A1,/,I3,/,I3)
0079 1000 FORMAT(' PLEASE SPECIFY DISK INPUT FILES WITH ONE '/'
*      ' ID CHARACTER <CR> (A1 FORMAT)',/)
*      ' THEN SPECIFY THE FIRST FILE NUMBER IN THE SEQ-'/

```

```

      *   ' UENCE <CR>, AND THE LAST FILE NUMBER <CR> (I3 FORMATS).'/)
1001  FORMAT(' INPUT WAVELENGTH SCALE (MICRONS/MICRON ON FILM),'/
      *   ' THEN <CR>, XSTEP (SECONDS OF ARC) BETWEEN SPECTRA, '/
      *   ' THEN <CR>, YUNIT (SECONDS OF ARC PER'/
      *   ' MICRON ON FILM, <CR>.(E10.0 FORMATS).'/)
1002  FORMAT(' INPUT WAVELENGTH OF FILM LEADING EDGE (MICRONS) <CR>,'/
      *   ' WAVELENGTH OF LINE CENTER <CR>.'/)
1003  FORMAT(' INPUT WAVELENGTHS OF ADDITIONAL LINES ON FILM'/
      *   ' THEN <CR>, OR <CR> IF NO MORE.'/)
1004  FORMAT(' INPUT POLARIZATION (RHC-511,LHC-256,LINEAR-ANGLE'/
      *   ' IN DEREES FROM RA 0-180) OF FIRST IN PAIR <CR>.'/
      *   ' OF SECOND IN PAIR <CR> (I3 FORMAT).'/)
1005  FORMAT(' INPUT IX (SLIT STEPS) FOR FIRST INPUT RECORD'/
      *   ' FROM BEGINNING OF THIS SCAN (>0 WITH FIRST SLIT CROSSING'/
      *   ' OF REGION SCAN = 1)(I3 FORMAT).'/)
1006  FORMAT(I4)
1007  FORMAT(' OPENING FILE ',10A1)
1008  FORMAT(' ERR:FORT0010',I2,' EXECUTION TERMINATED.')
1008  END

```

ROUTINES CALLED:

SETFIL, CURVE, READ, SPARM, PROCSS, DELETE

SWITCHES = /ON

BLOCK	LENGTH
MAIN	1394 (005344)*
BUFFER	2 (000004)

```

**COMPILER -- CORE**
      PHASE      USED  FREE
DECLARATIVES  00366 13425
EXECUTABLES   00967 12824
ASSEMBLY      01729 14979

```

```

0001      SUBROUTINE READ(NEWREC,N,ICODE,IBUF,END)
C      READ IS A RANDOM ACCESS READ STATEMENT FOR DENSITOMETER
C      DATA RECORDS. NEWREC SPECIFIES THE RECORD TO BE READ, AND
C      ICODE,N,IBUF(1) TO IBUF(N) ARE THE INFORMATION ON THAT RECORD. END
C      IS SET .TRUE. IF AN ATTEMPT HAS BEEN MADE TO READ BEYOND
C      THE LAST RECORD OF THE FILE.
0002      DIMENSION IBUF(1)
0003      LOGICAL END
0004      DATA IREC/0/
0005      END=.FALSE.
0006      IF(NEWREC=IREC) 30,25,10
0007      10 IREC=IREC+1
0008      DO 20 L1=IREC,NEWREC
0009      READ(1,END=60) ICODE,N,(IBUF(L2),L2=1,N)
0010      20 CONTINUE
0011      IREC=NEWREC
0012      25 RETURN
0013      30 REWIND 1
0014      IREC=0
0015      GO TO 10
0016      60 END=.TRUE.
0017      RETURN
0018      END

```

SWITCHES = /ON

BLOCK LENGTH
 READ 140 (000430)+

COMPILER --- CORE
 PHASE USED FREE
 DECLARATIVES 00366 13425
 EXECUTABLES 00527 13264
 ASSEMBLY 00996 15712


```

0001      SUBROUTINE SPARM(IBUF,ZUNIT,XSTEP,YUNIT,NPPL,NLPF,NFR,
      *      LPAT,ZCEN,ZSTEP)
C      SPARM EXAMINES THE PARAMETER RECORD OF THE INPUT DENSITOMETER
C      DATA AND RETURNS CERTAIN PARAMETERS ESSENTIAL FOR THE PROGRAM.
C      THESE PARAMETERS ARE:
C      NPPL=NUMBER OF POINTS PER LINE.
C      NLPF=NUMBER OF LINES PER FRAME=NUMBER OF YSTEPS IN SCAN.
C      LPAT=PATTERN=.TRUE. IF RASTER,.FALSE. IF EDGE.
C      NFR=NUMBER OF FRAMES=1 IF OLD DATA, 2 IF NEW DATA.
C      IN ADDITION SPARM SETS /FILM/ COMMON AND WRITES THE PARAMETER RECO
C      FOR THE OUTPUT TAPE.
0002      COMMON /FILM/SLTSTP(2),SMRSTP(2),SEFF,SLT(2),SMARG(2),ZEFF
0003      DIMENSION IBUF(1)
0004      DATA ICODE/2/,N/10/
0005      NPPL=IBUF(1)
0006      IDELX=IBUF(4)
0007      IDELY=IBUF(5)
0008      NLPF=IBUF(8)
0009      LPAT=IBUF(9)
0010      NFR=IBUF(12)/2+1
0011      IF(NFR.LT.1) NFR=1
0012      IF(NFR.GT.2) NFR=2
0013      DELX=IDELX
0014      DO 10 L1=1,2
0015          SLTSTP(L1)=SLT(L1)/DELX
0016          SMRSTP(L1)=SMARG(L1)/DELX
0017      10 CONTINUE
0018      YSTEP=YUNIT*FLOAT(IDELY)
0019      ZSTEP=ZUNIT*DELX
0020      SEFF=ZEFF/ZSTEP
0021      WRITE(8) ICODE,N,XSTEP,YPSTEP,ZSTEP,ZCEN,NLPF,NFR
0022      RETURN
0023      END

```

ROUTINES CALLED:

FLOAT

SWITCHES = /ON

BLOCK	LENGTH
SPARM	278 (001054)*
FILM	20 (000050)

```

**COMPILER --- CORE**
      PHASE      USED  FREE
DECLARATIVES 00366 13425
EXECUTABLES  00701 13090
ASSEMBLY     01141 15567

```

```

0001      SUBROUTINE PRCS$(INIT, IDAT, IPOL, IX0, IY, ZSTEP, NLNE, ZLNE,
      *      LPAT, NPPL, IFR, NFR)
C      PRCS$ EXAMINES DENSITOMETER DATA RECORDS AND PRODUCES
C      SEPERATED OUTPUT DATA RECORDS EACH CONTAINING INFORMATION
C      FOR A SINGLE SLIT CROSSING.
C      THIS IS A FOUR STEP PROCESS. FIRST, IF THE DATA IS A RASTER BACK
C      SCAN LINE THE DATA IS FLIPPED. SECOND, CR3DEG FINDS THE
C      EDGE OF THE FIRST SLIT CROSSING ON THE RECORD. THIRD, THE PROGRAM
C      WALKS THROUGH THE LINE SEPERATING OUT SLIT CROSSINGS AND
C      WRITES THESE ON A SCRATCH FILE. FOURTH, DATA IS WRITTEN ON TAPE
C      IN A CONVENIENT ORDER FOR LATER USE.
0002      DIMENSION IDAT(1), IPOL(1), ZLNE(1), ILNE(3), INTG(2)
0003      LOGICAL LPAT, LDIR, SKIP
0004      COMMON /FILM/SLT(2), SMARG(2), SEFF
0005      EQUIVALENCE (INTG(1), REALX)
0006      DATA ICODE/0/
0007      IF (INIT.EQ.-1) LDIR=.TRUE.
0008      IF (.NOT.LDIR) CALL FLIP(IDAT, NPPL)
0009      IF (INIT.EQ.-1) CALL CR3DEG(IDAT, ZERO, NPPL, NFR)
0010      BUF=ZERO-SMARG(IFR)
0011      L1=0
0012      20 L1=L1+1
0013      IF (NFR.NE.2) L2=L1-2+((L1-1)/2)
0014      IF (NFR.EQ.2) L2=IFR
0015      BUF=BUF+SMARG(L2)
0016      NBUF=BUF+.5+SLT(L2)
0017      IF (NBUF.GT.NPPL) GO TO 40
0018      IREC2=NFR*((L1-1)+IFR)
0019      IDAT(1)=SLT(L2)
0020      IF (SLT(L2).GT.SEFF) IDAT(1)=SEFF
0021      IDAT(2)=IX0+(L1-1)/2
0022      IDAT(3)=IY
0023      I1=SLT(L2)*SEFF
0024      IF (I1.LT.0) I1=0
0025      J1=BUF+1.
0026      J2=BUF+SLT(L2)
0027      REALX=ZER(IDAT, J1, J2, NLNE, ZLNE, ZSTEP, INIT)-FLOAT(I1/2)
0028      IDAT(4)=INTG(1)
0029      IDAT(5)=INTG(2)
0030      IDAT(6)=IPOL(L2)
0031      CALL PHOT(IDAT, J1, J2)
0032      I1=SLT(L2)
0033      L4=BUF
0034      DBUF=BUF-FLOAT(L4)
0035      IWAY=+1
0036      IF (DBUF.GE.0.) GO TO 23
0037      IWAY=-1
0038      DBUF=ABS(DBUF)
0039      23 DBUF I=1.-DBUF
0040      IF (SLT(L2).LE.SEFF) GO TO 22
0041      I1=SEFF
0042      L4=BUF+SEFF/2.
0043      22 J1=7
0044      J2=I1+6
0045      DO 30 L3=J1, J2

```

```

0046      L4=L4+1
0047      30 IDAT(L3)=DBUF*IDAT(L4+IWAY)+DBUF1*IDAT(L4)
0048      BUF=BUF+SLT(L2)
0049      WRITE(2,'IREC2) J2,(IDAT(L3),L3=1,J2)
0050      GO TO 20
0051      40 CONTINUE
0052      IF(IFR.LT.NFR) RETURN
0053      IF(IREC2-2*(IREC2/2).NE.0) IREC2=IREC2-1
0054      DO 50 L1=1,IREC2
0055      READ(2,'L1) I1,(IDAT(L2),L2=1,I1)
0056      WRITE(8) ICODE,I1,(IDAT(L2),L2=1,I1)
0057      50 CONTINUE
0058      IF(LPAT) LDIR=.NOT.LDIR
0059      RETURN
0060      END
    
```

ROUTINES CALLED:

FLIP , CR3DEG, ZER , FLOAT , PHOT , ABS

SWITCHES = /ON

BLOCK	LENGTH
PROSS 855	(003256)*
FILM 10	(000024)

COMPILER --*- CORE

PHASE	USED	FREE
DECLARATIVES	00366	13425
EXECUTABLES	00948	12843
ASSEMBLY	01513	15195

```
0001      SUBROUTINE FLIP(IDAT,N)
      C    FLIP TAKES AN ARRAY IDAT OF LENGTH N AND REVERSES THE ORDER OF
      C    THE ARRAY ELEMENTS.
0002      DIMENSION IDAT(1)
0003      NP1=N+1
0004      MID=N/2
0005      DO 10 L1=1,MID
0006      L1P=NP1-L1
0007      IHL0=IDAT(L1P)
0008      IDAT(L1P)=IDAT(L1)
0009      IDAT(L1)=IHL0
0010  10 CONTINUE
0011      RETURN
0012      END
```

SWITCHES = /ON

BLOCK	LENGTH
FLIP	107 (000326)*

COMPILER --- CORE

PHASE	USED	FREE
DECLARATIVES	00366	13425
EXECUTABLES	00527	13264
ASSEMBLY	00940	15768

```

-0001      FUNCTION ZER(IDAT,J1,J2,NLNE,ZLNE,ZSTEP,INIT)
      C     ZER DETERMINES THE DISTANCE FROM THE BEGINNING OF THE SLIT
      C     CROSSING TO LINE CENTER BY CALIBRATING OFF OF LINE EDGE AND ANY
      C     OTHER REFERENCE LINES IN THE FRAME.
0002      DIMENSION IDAT(1),ZLNE(1)
0003      LOGICAL INIT
0004      IF(NLNE.NE.0) GO TO 10
0005      IF(.NOT.INIT) RETURN
0006      ZER=(ZLNE(2)-ZLNE(1))/ZSTEP
0007      RETURN
0008  10 CONTINUE
0009      ZC=1.
0010      DO 40 L1=1,NLNE
0011      ILNE=(ZLNE(L1+2)-ZLNE(1))/ZUNIT-10.
0012      J2P=ILNE+50
0013      J2M=ILNE+10
0014      IDIFM=0
0015      DO 20 L2=ILNE,J2P
0016      IDIF=IDAT(L2+J1+1)-IDAT(L2+J1)
0017      IF(IDIF.LT.IDIFM) GO TO 20
0018      IDIFM=IDIF
0019      L2M=L2
0020  20 CONTINUE
0021  40 ZC=ZC+FLOAT(L2M)*ZSTEP/(ZLNE(L1+2)-ZLNE(1))
0022      ZER=ZC/FLOAT(NLNE+1)*(ZLNE(2)-ZLNE(1))/ZSTEP
0023      RETURN
0024      END

```

ROUTINES CALLED:

FLOAT

SWITCHES = /ON

BLOCK	LENGTH
ZER	326 (001214)*

```

**COMPILER --*- CORE**
      PHASE      USED  FREE
DECLARATIVES 00446 13345
EXECUTABLES  00607 13184
ASSEMBLY     01132 15576

```

```

0001      SUBROUTINE CR3DEG(IDAT,ZER,NPPL,NFR)
C      CR3DEG FINDS THE EDGE OF THE FIRST SLIT CROSSING IN THE DATA A
C      IDAT. IT LOOKS FOR A SLOPE IN FILM DENSITIES THAT IS GREATER
C      THAN ICONTR OVER TEN POINTS.
0002      DIMENSION IDAT(1)
0003      COMMON /DENS/ICONTR
0004      IDIFF(N,M)=(IDAT(N)-IDAT(N+M))*ISGN
0005      ISGN=1
0006      IF(NFR.EQ.2) ISGN=-1
0007      LIMITN=NPPL/2
0008      DO 20 L1=1,LIMITN
0009      IDIF=IDIFF(L1,10)
0010      IF(IDIF.GT.ICONTR) GO TO 30
0011  20 CONTINUE
0012      WRITE(6,'1000')
0013      STOP 2
0014  30 NZER=L1+1
0015      IDIFM=0
0016      LIM=0
0017      DO 40 L1=1,10
0018      IDIF=IDIFF(NZER+L1,1)
0019      IF(IDIF.LT.IDIFM) GO TO 40
0020      IDIFM=IDIF
0021      LIM=L1
0022  40 CONTINUE
0023      NZER=NZER+L1M
0024      ZER=NZER
0025      RETURN
0026  1000 FORMAT(' ERR:CR3DEG COULD NOT FIND EDGE. EXECUTION TERMINATED')
0027      END

```

SWITCHES = /ON

BLOCK	LENGTH
CR3DEG 286	(001074)*
DENS 1	(000002)

```

**COMPILER --- CORE**
  PHASE      USED  FREE
DECLARATIVES 00366 13425
EXECUTABLES  00616 13175
ASSEMBLY     01245 15463

```

```

0001      SUBROUTINE CURVE
C          PHOTOMETRY FINDS 1024-EXPOSURE AS FCN OF TRANS./2
C          USES HARVARD PHOTOMETRY CURVE OF JUNE, 1971 (5 PERCENT)
C          T. TARBELL NOV. 19, 1973
0002      COMMON/CCURV/IPHOT(513),KSMOOT
0003      DATA C1,C2,C3,A0,A1,A2,SLOPE,EINT,B0,B1,B2,S2,E2/
1          2.017,0.408,0.1115,101.61,218.26,-13.1,161.63,162.53,
2          44.46,-7740.5,69480.0,1571.07,-123.87/
0004      DO 10 I=1,513
0005      F=FLOAT(I)
0006      F=ALOG(500./F)+0.1
0007      IF (F.GT.C1) GO TO 1
0008      IF (F.GT.C2) GO TO 2
0009      IF (F.GT.C3) GO TO 3
0010      IPHOT(I) = S2*F + E2
0011      IF(IPHOT(I).LT.1) IPHOT(I)=1
0012      GO TO 9
0013      1 IPHOT(I) = SLOPE*F + EINT
0014      IF(IPHOT(I).GT.1023) IPHOT(I) = 1023
0015      GO TO 9
0016      2 IPHOT(I) = A0 + A1*F + A2*F*F
0017      GO TO 9
0018      3 IPHOT(I) = B0 + SQRT(B1 + B2*F)
0019      9 IPHOT(I) = 1.40*IPHOT(I)
0020      IPHOT(I) = MAX0(1,1024-IPHOT(I))
0021      IPHOT(I)=MIN0(1023,IPHOT(I))
0022      10 CONTINUE
0023      WRITE (6,100)
0024      100 FORMAT(/' TYPE 1 FOR SMOOTHING, 0 FOR NONE'/)
0025      READ (6,101) KSMOOT
0026      101 FORMAT (I2)
0027      RETURN
0028      END

```

ROUTINES CALLED:

FLOAT , ALOG , SQRT , MAX0 , MIN0

SWITCHES = /ON

BLOCK	LENGTH
CURVE 422	(001514)*
CCURV 514	(002004)

COMPILER --*-+ CORE

PHASE	USED	FREE
DECLARATIVES	00537	13254
EXECUTABLES	00687	13104
ASSEMBLY	01397	15311

```

0001      SUBROUTINE PHOT(IDAT,J1,J2)
          C      CONVERTS FROM TRANSMISSION TO (1024*EXPOSURE)
          C      T. TARRELL NOV. 19, 1973
0002      DIMENSION IDAT(1), RAT(10)
0003      COMMON /WEDGE/NSTEPS,NOWEG,ITRANS(11),INT(11)
0004      COMMON /CCURV/IPHOT(513),KSMOOT
0005      LOGICAL NOWEG
0006      DO 100 L=J1,J2
0007      LL=IDAT(L)/2
0008      IF(LL.LE.0) LL=1
0009      IF(LL.GE.513) LL=513
0010      100 IDAT(L)=IPHOT(LL)
0011      IF(KSMOOT.EQ.0) RETURN
0012      M2=J1+2
0013      M3=J2-2
0014      DO 200 L1=M2,M3
0015      200 IDAT(L1)=0.4*(IDAT(L1)+0.5*(IDAT(L1+1)+IDAT(L1-1)+
          1      0.5*(IDAT(L1+2)+IDAT(L1-2))))
0016      RETURN
0017      END

```

SWITCHES = /ON

BLOCK	LENGTH
PHOT	241 (000742)*
WEDGE	24 (000060)
CCURV	514 (002004)

```

**COMPILER ---- CORE**
  PHASE      USED  FREE
DECLARATIVES 00540 13242
EXECUTABLES  00607 13184
ASSEMBLY     01098 15610

```


0001

BLOCK DATA

C ALL THE PROGRAM CONTROL PARAMETERS ARE SET HERE. THESE PARA-
C ARE DESCRIBED:
C /FILM/=CONTAINS INFORMATION ON FILM DIMENSIONS.
C SLT(1 AND 2)=MICRONS ACROSS SLIT FOR FIRST (1) AND SECOND (2)
C ORTHOGONAL POLARIZATIONS.
C SMARG(1 AND 2)=MICRONS IN MARGIN PRIOR TO SLIT CROSSING.
C ZEFF=NUMBER OF MICRONS OF USEFUL DATA AROUND LINE CENTER.
C /DENS/=ICONTR=MINIMUM CONTRAST BETWEEN SLIT EDGE AND MARGIN.
C /WEDGE/=STEP WEDGE DATA.
C NOWEG=.TRUE. IF NO WEDGE INFORMATION IS AVAILABLE
C NSTEPS=NUMBER OF WEDGE STEPS.
C ITRANS(K)=TRANSMISSION OF K-TH STEP IN WEDGE.
C INT(K)=CORRECT TRANSMISSION OF K-TH STEP IN WEDGE.

0002

COMMON /FILM/DUM1(5),SLT(2),SMARG(2),ZEFF

* /DENS/ICONTR
* /BUFFER/IC,IBUF(10240)
* /WEDGE/NSTEPS,NOWEG,ITRANS(11),INT(11)

0003

LOGICAL NOWEG

0004

DATA ICONTR/200/

0005

DATA SMARG,SLT/3527.,685.,2644.,2644./,ZEFF/0.8/

0006

DATA NSTEPS/11/,ITRANS/0.100,200,300,400,500,600,

* 700,800,900,1000/,INT/0.100,200,300,400,500,600,700,
* 800,900,1000/

0007

DATA NOWEG/.TRUE./

0008

END

SWITCHES = /ON

BLOCK	LENGTH
DATA	0 (000000) *
FILM	20 (000050)
DENS	1 (000002)
BUFFER	10241 (050002)
WEDGE	24 (000060)

COMPILER --*- CORE
PHASE USED FREE
DECLARATIVES 00366 13425
EXECUTABLES 00570 13221
ASSEMBLY 00860 15848

```

C PASS 3 DATA FROM TAPE TO FINAL TAPE FORM
C THIS PROGRAM TAKES A SEQUENCE OF INPUT FILES LABELED <CHAR>1.DA
C THROUGH <CHAR>M.DAT OFF OF TAPE WHICH HAS BEEN MADE BY PASS2 AND
C DETERMINES THE MAGNETIC FIELDS B, FIELD ANGLES (GAM,PHI), AND
C VELOCITIES V FOR ALL THE POINTS IN THE REGION. PASS3 STORES THE
C RESULTS OF ITS ANALYSIS IN A SCRATCH FILE UNTIL COMPLETION OF
C THE ANALYSIS. AT COMPLETION OF THE ANALYSIS PASS3 REQUESTS AN OUT
C TAPE AND AFTER CONTINUATION PUTS THE DATA ON TAPE. THE FINAL TAPE
C RECORDS BEGIN WITH TWO KEY WORDS ICODE AND N. ICODE
C DEFINES THE NATURE OF THE RECORD AND N THE LENGTH OF THE BUFFER TO
C FOLLOW. THE OUTPUT RECORDS ARE CODED AS FOLLOWS:
C ICODE=1,N,IBUF=LABEL OF LENGTH N CHARACTERS.
C ICODE=2,N,IBUF(1&2)=XSTEP=SECONDS OF ARC BETWEEN X STEPS-REAL).
C IBUF(3&4)=YSTEP=SECONDS OF ARC BETWEEN Y STEPS-REAL).
C IBUF(5)=NX=NUMBER OF FIELD UNITS IN X DIMENSION.
C IBUF(6)=NY=NUMBER OF FIELD UNITS IN Y DIMENSION.
C ICODE=0,N,IBUF(1)=IX=X STEP POSITON OF THIS DATA POINT.
C IBUF(2)=IY=Y STEP POSITON OF THIS DATA POINT.
C IBUF(3&4)=V=VELOCITY AT THIS POINT(KM/SEC, +=DOWN)-(REAL).
C IBUF(5&6)=B=FIELD STRENGTH AT THIS POINT(GAUSS)-(REAL).
C IBUF(7&8)=GAM=FIELD ANGLE FROM LINE OF SIGHT(DEGREES)-(REAL).
C IBUF(9&10)=PHI=FIELD ANGLE FROM X SCAN DIRECTION(DEGREES)-(REAL).
C LARRY NOVEMBER - 21-AUG-73

```

```

0001 COMMON /BUFFER/IBUFC(512)
0002 DIMENSION IBUF(256,2),IBUF1(256),IBUF2(256),IEXT(5),LABEL(32)
0003 LOGICAL*1 CHAR(2),IFIL(10),IFILNM(4)
0004 EQUIVALENCE (IBUFC(1),IBUF(1,1),IBUF1(1)),(IBUFC(257),IBUF2(1)),
* (IFIL(1),IEXT(1)),(IFIL(3),IFILNM(1)),(CHAR(1),ICHAR)
0005 DATA IEXT/'00','00','00','.D','AT'/,ICHAR/' 0'/
0006 CALL SETFIL(2,'SCRT(H.TMP',IERR,'DK',0,0,'233,2)
0007 WRITE(6,1000)
0008 READ(6,901) NX,NY
0009 WRITE(6,1001)
0010 READ(6,900) IFIL(1)
0011 L1=0
0012 1 L1=L1+1
0013 ENCODE(4,1003,IFILNM) L1
0014 DO 3 L2=1,4
0015 3 IF(IFILNM(L2).EQ.CHAR(1)) IFILNM(L2)=CHAR(2)
0016 CALL SETFIL(8,IFIL)
0017 5 CONTINUE
0018 READ(8,END=100,ERR=200) ICODE,N,(IBUF1(L2),L2=1,N)
0019 I1=ICODE+1
0020 GO TO (10,20,40),I1
0021 10 READ(8,END=100) ICODE,N,(IBUF2(L2),L2=1,N)
0022 CALL PROCSS(IBUF,NX,ZSTEP)
0023 GO TO 5
0024 20 IF(L1.NE.1) GO TO 5
0025 DO 30 L2=1,N
0026 30 LABEL(L2)=IBUF1(L2)
0027 LABEL(32)=N
0028 GO TO 5
0029 40 IF(L1.NE.1) GO TO 5
0030 CALL SPARM3(IBUF1,NY,XSTEP,YSTEP,ZSTEP)
0031 CALL IO1(3,IBUF1,NX,NY,IERR)

```

```

0032      NX=0
0033      GO TO 5
0034      100 END FILE 8
0035      GO TO 1
0036      200 CONTINUE
0037      DO 210 L1=1,10
0038      CALL BELL
0039      DO 210 L2=1,10000
0040      210 CONTINUE
0041      WRITE(6,1002)
0042      PAUSE 1
0043      ICODE=1
0044      N=LABEL(32)
0045      WRITE(8) ICODE,N,(LABEL(L1),L1=1,N)
0046      ICODE=2
0047      N=6
0048      WRITE(8) ICODE,N,XSTEP,YSTEP,NX,NY
0049      ICODE=0
0050      N=10
0051      DO 250 L1=1,NY
0052      DO 250 L2=1,NX
0053      CALL IO1(1,IBUF1,L2,L1,IREF)
0054      V=FLOAT(IBUF1(IREF+1))/256.
0055      B=FLOAT(IBUF1(IREF+2))/4.
0056      GAM=FLOAT(IBUF1(IREF+3))/128.
0057      PHI=FLOAT(IBUF1(IREF+4))/128.
0058      250 WRITE(8) ICODE,N,L2,L1,V,B,GAM,PHI
0059      END FILE 8
0060      END FILE 2
0061      CALL DELETE
0062      900 FORMAT(A1)
0063      901 FORMAT(I3)
0064      1000 FORMAT('! SPECIFY THE MAXIMUM NUMBER OF PICTURE!',
*      ' ELEMENTS IN X <CR>.'/
*      ' AND THE MAXIMUM NUMBER OF PICTURE.',
*      ' ELEMENTS IN Y <CR>.'/ (I3 FORMATS).'/)
0065      1001 FORMAT('! SPECIFY ONE LETTER WHICH IDENTIFIES TAPE!',
*      ' INPUT FILES<CR>.'/)
0066      1002 FORMAT('! MOUNT TAPE FOR OUTPUT OF MAGNETOGRAM AND ASSIGN.'/
*      ' TAPE FILE TO UNIT 8, THEN CONTINUE.'/)
0067      1003 FORMAT(I4)
0068      END

```

ROUTINES CALLED:

SETFIL, PRCS5, SPARM3, IO1, BELL, FLOAT, DELETE

SWITCHES = /ON

BLOCK	LENGTH
MAIN	995 (003706)*
BUFFER	512 (002000)

COMPILER ---- CORE

PHASE	USED	FREE
DECLARATIVES	00366	13425
EXECUTABLES	00911	12880
ASSEMBLY	01629	15079

```

0001      SUBROUTINE SPARM3(BUF,NY,XSTEP,YSTEP,ZSTEP)
C      SPARM3 INTERPRETS THE PARAMETER RECORD ON THE INPUT
C      TAPE AND RETURNS TO THE PROGRAM THE FOLLOWING PARAMETERS:
C      XSTEP=SECONDS OF ARC PER X INCREMENT.
C      YSTEP=SECONDS OF ARC PER Y INCREMENT.
C      ZSTEP=MICRONS WAVELENGTH BETWEEN DATA POINTS.
C      NY=TOTAL NUMBER OF STEPS IN Y DIMENSION.
C      DVEL=KM/SEC PER MICRON WAVELENGTH OF RED SHIFT.
C      DFLD=GAUSS PER MICRON WAVELENGTH OF COMPONENT SHIFT
0002      DIMENSION BUF(4),INTG(2)
0003      COMMON /LNPRM/DVEL,DFLD,GFAC
0004      EQUIVALENCE (INTG(1),REALX)
0005      XSTEP=BUF(1)
0006      YSTEP=BUF(2)
0007      ZSTEP=BUF(3)
0008      ZCEN=BUF(4)
0009      REALX=BUF(5)
0010      NY=INTG(1)
0011      DVEL=2.998E5/ZCEN
0012      DFLD=2.142E10/GFAC/ZCEN**2
0013      RETURN
0014      END

```

SWITCHES = /ON

BLOCK	LENGTH
SPARM3 135	(000416)*
LNPRM 6	(000014)

```

**COMPILER --- CORE**
  PHASE      USED  FREE
DECLARATIVES 00366 13425
EXECUTABLES  00628 13163
ASSEMBLY     00985 15723

```

```

0001      SUBROUTINE PRCSS(IBUF,NX,ZSTEP)
C      PRCSS DETERMINES THE FIELD STRENGTH, DIRECTION, AND VELOCITY
C      FROM THE RAW DATA IBUF USING THE THREE ROUTINES PARM,FIELD, AND
C      PRMFFT.
0002      DIMENSION PK(2,15),IBUF(256,2),IZ0(2),DAT(256),IPOL(2),
*      PKH(15),DATA(256,2)
0003      EQUIVALENCE (IZ0(1),Z0)
0004      COMMON /FILTER/NSAMP /FFT/NOFFT,BMIN,GBND
0005      LOGICAL NOFFT
0006      DATA NTIME/-1/,NANG/-1/
0007      DO 50 L1=1,2
0008      N=IBUF(1,L1)
0009      IX=IBUF(2,L1)
0010      IY=IBUF(3,L1)
0011      IZ0(1)=IBUF(4,L1)
0012      IZ0(2)=IBUF(5,L1)
0013      IPOL(1)=IBUF(6,L1)
0014      DO 10 L2=1,N
0015      10 DAT(L2)=IBUF(L2+6,L1)
0016      DO 30 L2=1,N
0017      30 DATA(L2,L1)=DAT(L2)
0018      50 CONTINUE
0019      IF((IX.EQ.1).AND.(IY.EQ.1)) NTIME=NTIME+1
0020      IF((IX.EQ.1).AND.(IY.EQ.1).AND.(IPOL(1).LT.256)) NANG=NANG+1
0021      IF((NTIME.LT.0).OR.((IPOL(1).LT.256).AND.(NANG.LT.0))) RETURN
0022      IF(IX.GT.NX) NX=IX
0023      CALL IO1(1,IBUF,IX,IY,IREC)
0024      GAM=FLOAT(IBUF(IREC+3,1))/128.
0025      PHI=0.
0026      CALL PRMFFT(DATA,N,ZSTEP,Z0,V,B,GAM,CONIN)
0027      B=INT(B/10.)
0028      B=B*10.+CONIN/100.
0029      I1=V*256.
0030      I2=B*4.
0031      I3=GAM*128.
0032      I4=PHI*128.
0033      IBUF(IREC+1,1)=(NTIME*IBUF(IREC+1,1)+I1)/(NTIME+1)
0034      IF(IPOL(1).GE.256) IBUF(IREC+2,1)=I2
0035      65 IF(IPOL(1).GE.256) IBUF(IREC+3,1)=I3
0036      IF(IPOL(1).LT.256) IBUF(IREC+4,1)=
*      (NANG*IBUF(IREC+4,1)+I4)/(NANG+1)
0037      CALL IO1(2,IBUF,IX,IY,IREC)
0038      RETURN
0039      END

```

ROUTINES CALLED:

IO1 , FLOAT , PRMFFT , INT

SWITCHES = /ON

BLOCK	LENGTH
PRCSS	2244 (010610)*
FILTER	1 (000002)
FFT	5 (000012)

COMPILER --- CORE
 PHASE USED FREE
 DECLARATIVES 00637 13154
 EXECUTABLES 00847 12944

```

0001      SUBROUTINE IO1(MODE,IBUF,IX,IY,IREC)
C      IO SCRATCH IS HANDLED BY IO1.  THREE MODES ARE POSSIBLE, MODE=
C      READS INTO IBUF THE REQUESTED COORDINATE IX,IY DATA.  THAT DATA
C      IS CONTAINED IN THE FOUR RECORDS FOLLOWING THE IREC ELEMENT IN ARH
C      IBUF.  MODE=2 WRITES ON SCRATCH IBUF.  MODE=3 INITIALIZES THE FILE
C      TO ACCEPT IX BY IY ELEMENTS.
0002      DIMENSION IBUF(1)
0003      X=IX
0004      Y=IY
0005      GO TO (100,200,300),MODE
0006      100 CONTINUE
0007      IF((IX.LE.0).OR.(IY.LE.0)) RETURN
0008      REC=(X-1.)*YDIM+Y
0009      ISEG=(REC-1.)/64.+1.
0010      IREC=(REC-FLOAT(ISEG-1)+64.+1.)*4.
0011      READ(2,ISEG) (IBUF(L1),L1=1,256)
0012      RETURN
0013      200 CONTINUE
0014      IF((IX.LE.0).OR.(IY.LE.0)) RETURN
0015      REC=(X-1.)*YDIM+Y
0016      ISEG=(REC-1.)/64.+1.
0017      WRITE(2,ISEG) (IBUF(L1),L1=1,256)
0018      RETURN
0019      300 CONTINUE
0020      YDIM=Y
0021      REC=X+Y
0022      ISEG=(REC-1.)/64.+1.
0023      DEFINE FILE 2(ISEG,256,U,IVAR)
0024      WRITE(2,'1,ERR=350) IDUMMY
0025      RETURN
0026      350 END FILE 2
0027      CALL DELETE
0028      WRITE(6,1000) IREC
0029      STOP 1
0030      1000 FORMAT('1,ERR=350',I2,' EXECUTION TERMINATED.1')
0031      END

```

ROUTINES CALLED:
 FLOAT , DELETE

SWITCHES = /ON

BLOCK	LENGTH
IO1	383 (001376)*

COMPILER --- CORE

PHASE	USED	FREE
DECLARATIVES	00366	13425
EXECUTABLES	00607	13184
ASSEMBLY	01300	15408

```

SUBROUTINE PRMFFT(DATA,N,ZSTEP,Z0,V,B,GAMF,CONIN)
MASTER CALLING ROUTINE FOR FT ANALYSIS OF ZEEMAN PROFS
NEEDS CONCEN,ZER,POLYAN
C
C
0001  DIMENSION DATA(512)
0002  DIMENSION GOLD(2,9)
0003  COMMON /LNPRM/DVEL,DFLD,GFAC
0004  COMMON /RESULTS/X1,COSGAM,GAM,SLOPE,IZERO,NORAT,
0005  1  PLCOF(5),ACOF,BCOF,CCOF,ALOR
0006  COMMON /STUFF/CONT(4),NCONT(4),COMT(2),ISTART(4),
1  IEND(4),SUM(140),PHD,DELNC,SHIFT,NSHIFT,SCHIFF,
2  NCENT,LCENT,SCALE,P,Q,LIMIT,AREA
0007  DATA DOPVEL/0./
0008  107  FORMAT (' N, P, Q ',I4,2F10.4)
0009  DO 9 I=1,N
0010  J=I+256
0011  DATA(J)=1024.-DATA(J)
0012  DATA(I)=1024.-DATA(I)
0013  9  CONTINUE
0014  101  FORMAT(' ',10F10.0)
0015  M3=256+N
C  PRINT 101, (DATA(I),I=1,N)
C  PRINT 101, (DATA(I),I=257,M3)
0016  DO 7 I=1,2
0017  DO 7 J=1,9
0018  7  GOLD(I,J)=0.
0019  V=0.
0020  LCENT=0
0021  NSHIFT=256
0022  NCENT=Z0+0.5+DOPVEL
0023  DOPVEL=0.
0024  P=0.05*(DATA(7)-DATA(NCENT))
0025  Q=1.1*P
0026  ISTART(1)=Z0/2+1
0027  ISTART(2)=(N+Z0)/2-1
0028  ISTART(3)=ISTART(1)+256
0029  ISTART(4)=ISTART(2)+256
0030  IEND(1)=6
0031  IEND(2)=N-5
0032  IEND(3)=262
0033  IEND(4)=251+N
0034  CALL CONTNA(DATA,N)
0035  5  CONTINUE
0036  COMT(1)=AMAX1(CONT(1),CONT(2))
0037  CONIN=COMT(1)
0038  COMT(2)=AMAX1(CONT(3),CONT(4))
C  ****FUDGE TO THROW OUT UNDEREXPOSED FRAMES ****
0039  DEXP=110.
0040  IF(ABS(COMT(2)-DATA(NCENT+NSHIFT)).LE.DEXP)
1  GO TO 1054
0041  M=0
0042  DO 21 J=1,257,256
0043  M=M+1
0044  M3=J+N-1
0045  DO 21 L=J,M3
0046  21  DATA(L)=COMT(M)-DATA(L)

```

```

0047      DMAX1=0.
0048      M5=NCONT(1)
0049      M15=NCONT(2)
0050      DO 2300 I=M5,M15
0051      IF(DATA(I).LE.DMAX1) GO TO 2300
0052      MM1=I
0053      DMAX1=DATA(I)
0054 2300    CONTINUE
0055      DMAX1=DSUM(MM1-2,MM1+2,DATA)
0056      DMAX2=0.
0057      M5=NCONT(3)
0058      M15=NCONT(4)
0059      DO 2301 I=M5,M15
0060      IF(DATA(I).LE.DMAX2) GO TO 2301
0061      MM1=I
0062      DMAX2=DATA(I)
0063 2301    CONTINUE
0064      DMAX2=DSUM(MM1-2,MM1+2,DATA)
0065      P=ALOG(DMAX2)/ALOG(DMAX1)
0066      IF(ABS(1.-P).GE.0.10) GO TO 2305
0067      DO 2302 I=1,256
0068 2302    DATA(I)=ABS(DATA(I))*P
0069      GO TO 2307
0070 2305    PRINT 2306,P
0071 2306    FORMAT (' P WRONG 'F10.2)
0072 2307    CONTINUE
0073      CALL COMBIN(DATA,N)
0074      DOPVEL=LCENT+PHD-20
C *** FUDGE FACTOR FOR 32 MICRON DATA PTS..5250
0075      V=0.2528*DOPVEL
0076      LA=LCENT-NCONT(1)
0077      LB=IABS(NCONT(2)-LCENT)
0078      LC=LCENT+NSHIFT-NCONT(3)
0079      LD=IABS(NCONT(4)-LCENT-NSHIFT)
0080      LIM=MIN0(LA,LB,69)-1
0081      DO 10 KMODE=1,2
0082      LIMIT=LIM-1
0083      NORIG=LIMIT+1
0084 100    FORMAT (' ',7I6)
0085      DO 964 L=1,LIMIT
0086      SUM(NORIG+L)=DATA(LCENT+L)
0087 964    SUM(NORIG-L)=DATA(LCENT-L)
0088      SUM(NORIG)=DATA(LCENT)
0089      CALL POLYAN(SUM,140,NORIG,LIMIT,PHD,KMODE)
C *** FUDGE FACTOR FOR 32 MICRON INTERVAL 5250 ***
0090      GOLD(KMODE,4)=SLOPE*113.02
0091      IF(IZERO.EQ.0) GO TO 71
0092      GOLD(KMODE,1)=X1
0093      GOLD(KMODE,2)=X1+113.02
0094      IF(NORAT.NE.0) GO TO 71
0095      GOLD(KMODE,3)=COSGAM
0096      GOLD(KMODE,5)=COSGAM*GOLD(KMODE,2)
0097      GOLD(KMODE,6)=PLCOF(2)
0098      GOLD(KMODE,7)=(ABS(PLCOF(3))+ABS(PLCOF(4))+ABS(PLCOF(5)))
0099      1 *0.33333/ABS(PLCOF(2))
      IF(ABS(COSGAM).LT.0.7) GO TO 71

```



```

0100      CGM=GOLD(KMODE,4)/GOLD(KMODE,2)
0101      GOLD(KMODE,3)=0.33*COSGAM+0.67*CGM
0102      GOLD(KMODE,5)=GOLD(KMODE,3)+GOLD(KMODE,2)
0103      71      CONTINUE
0104      LCENT=LCENT+NSHIFT
0105      LIM=MIN0(LC,LD,69)-1
0106      10      CONTINUE
0107      LCENT=LCENT-2*NSHIFT
C          PRINT 100,NCONT,CONT,LCENT,DOPVEL,V,
C          1      ((GOLD(I,J),J=1,8),I=1,2)
0108      1054     CONTINUE
0109      CALL BGOOD(GOLD,B,GAMF)
0110      100      FORMAT('0,4I5,4F10.1/, LCENT, DOPVEL, V '
0111      1      I4,2F10.3/(' '3F10.1,5F10.3)//)
0111      RETURN
0112      END

```

ROUTINES CALLED:

CONTNA, AMAX1, ABS, DSUM, ALOG, COMBIN, IABS
 MIN0, POLYAN, BGOOD

SWITCHES = /ON

BLOCK	LENGTH
PRMFFT	1532 (005770)*
LNPRM	6 (000014)
RESLTS	28 (000070)
STUFF	324 (001210)

COMPILER --- CORE

PHASE	USED	FREE
DECLARATIVES	00823	12968
EXECUTABLES	01247	12544
ASSEMBLY	02067	14641

```

0001      SUBROUTINE CONTNA(DATA,NPTS)
C      FINDS CONTINUA AT NCONT(I),CONT INT = CONT(I)
C      ISTART(4) CONTAINS INDICES OF FIRST PT. TESTED. IEND(4), LAST PT.
C      P IS MAX. ALLOWED EXPOSURE CHANGE AT +-2 PTS. FROM PY. TESTED.
C      Q, LIKEWISE AT +-4 PTS. AWAY.
C      MODIFIED 1/25/74 CONT(I) FOUND BY LOOKING FOR FLATNESS
C      NCONT(I) FROM MAX. SIGNAL NEAR FRAME EDGE
C      NEEDS FCN. DSUM
C      T. TARBELL JULY 12,1971.

0002      COMMON /STUFF/CONT(4),NCONT(4),COMT(2),ISTART(4),
1      IEND(4),SUM(70),DIT(70),PHD,DELNC,SHIFT,NSHIFT,SCHIFF,
2      NCENT,LCENT,SCALE,P,Q,LIMIT,AREA
0003      DIMENSION DATA(512)
0004      DIF(M,N)=ABS(DATA(M)-DATA(N))
0005      DO 10 I=1,4
0006      J=ISTART(I)
0007      K=1
0008      IF((I.EQ.1).OR.(I.EQ.3)) K=-1
0009      CTEMP=100.
0010      NTEMP=0
0011      1 J=J+K
0012      IF(((I.EQ.1).OR.(I.EQ.3)).AND.(J.LE.IEND(I))) GO TO 11
0013      IF(((I.EQ.2).OR.(I.EQ.4)).AND.(J.GE.IEND(I))) GO TO 11
0014      A=DIF(J,2)
0015      IF(A.GT.P) GO TO 1
0016      B=DIF(J,-2)
0017      IF(B.GT.P) GO TO 1
0018      C=DIF(J,4)
0019      IF(C.GT.Q) GO TO 2
0020      D=DIF(J,-4)
0021      IF(D.GT.Q) GO TO 2
0022      NCONT(I)=J
0023      CONT(I)=DSUM(J-4,J+4,DATA)/9.
0024      GO TO 10
0025      2 C=AMAX1(A,B)
0026      IF (NTEMP.EQ.0) GO TO 3
0027      IF ((C-CTEMP).GT.5.) GO TO 1
0028      IF (DATA(J).LT.DATA(NTEMP)) GO TO 1
0029      3 CTEMP=C
0030      NTEMP = J
0031      GO TO 1
0032      11 IF (NTEMP.GT.0) GO TO 12
0033      NCONT(I)=0
0034      CONT(I)=0.0
0035      GO TO 10
0036      12 NCONT(I)=NTEMP
0037      CONT(I)=DSUM(NTEMP-2,NTEMP+2,DATA)/5.
0038      10 CONTINUE
0039      L1=NCENT/4+1
0040      L2=3*(NPTS-NCENT)/4+NCENT-1
0041      NP=NPTS-1
0042      M3=NSHIFT+1
0043      DO 30 J=1,M3,NSHIFT
0044      L1=L1+J-1

```

```

0045      L2=L2+J-1
0046      NP=NP+J-1
0047      K1=J+2
0048      K2=NP
0049      T1=0.
0050      T2=0.
0051      20      IF(DATA(K1).LT.T1) GO TO 22
0052      N1=K1
0053      T1=DATA(K1)
0054      22      IF(DATA(K2).LT.T2) GO TO 23
0055      N2=K2
0056      T2=DATA(K2)
0057      23      K1=K1+1
0058      K2=K2+1
0059      IF((K1.GT.L1).AND.(K2.LT.L2)) GO TO 24
0060      IF(K1.GT.L1) K1=L1
0061      IF(K2.LT.L2) K2=L2
0062      GO TO 20
0063      24      IF(J.GT.1) GO TO 25
0064      NCONT(1)=N1
0065      NCONT(2)=N2
0066      IF((CONT(1).NE.0.).OR.(CONT(2).NE.0.)) GO TO 30
0067      CONT(1)=0.333*(DATA(N1)+DATA(N1-1)+DATA(N1+1))
0068      CONT(2)=0.333*(DATA(N2)+DATA(N2-1)+DATA(N2+1))
0069      GO TO 30
0070      25      NCONT(3)=N1
0071      NCONT(4)=N2
0072      IF((CONT(3).NE.0.).OR.(CONT(4).NE.0.)) GO TO 30
0073      CONT(3)=0.333*(DATA(N1)+DATA(N1-1)+DATA(N1+1))
0074      CONT(4)=0.333*(DATA(N2)+DATA(N2-1)+DATA(N2+1))
0075      30      CONTINUE
0076      RETURN
0077      END

```

ROUTINES CALLED:

ABS , DSUM , AMAX1

SWITCHES = /ON

BLOCK	LENGTH
CONTNA 1009	(003742)*
STUFF 324	(001210)

COMPILER --- CORE

PHASE	USED	FREE
DECLARATIVES	00714	13077
EXECUTABLES	00927	12864
ASSEMBLY	01545	15163

```

0001      SUBROUTINE COMBIN(DATA,NPTS)
          C      FINDS CENTRAL WAVELENGTH OF PROFILE USING ZEROth AND
          C      FIRST MOMENTS OF BOTH POLARIZATIONS
          C      CENTER IS AT K+PHD
          C      TDT 1/17/74
0002      DIMENSION DATA(512)
0003      EQUIVALENCE (K,LCENT)
0004      COMMON/STUFF/CONT(4),NCONT(4),COMT(2),ISTART(4),
          1      IEND(4),SUM(70),DIF(70),PHD,DELNC,SHIFT,NSHIFT,
          2      SCHIFF,NCENT,LCENT,SCALE,P,Q,LIMIT,AREA
0005      ST=0.
0006      NIT=0
0007      NPT=NPTS+NSHIFT
0008      PHD=0.
0009      1      NIT=NIT+1
0010      LA=NCENT-NCONT(1)
0011      LB=IABS(NCENT-NCONT(2))
0012      LC=IABS(NCENT-NPTS+1)
0013      LIM=MIN0(LA,LB,LC)-2
0014      W=0.1/FLOAT(LIM)
0015      AM0=COSTR(DATA,512,NCENT,LIM,W,ST,PHD,1)
0016      AM1=ST/W
0017      NCEN=NCENT+NSHIFT
0018      LA=NCEN-NCONT(3)
0019      LB=IABS(NCEN-NCONT(4))
0020      LIM=MIN0(LA,LB,LC)-2
0021      W=0.1/FLOAT(LIM)
0022      BM0=COSTR(DATA,512,NCEN,LIM,W,ST,PHD,1)
0023      IF(ABS(AM0-BM0).GT.3000.) GO TO 11
0024      BM1=ST/W
0025      DELNC=0.5*(BM1/BM0+AM1/AM0)
0026      CENA=NCENT+PHD
0027      CEN=CENA+DELNC
0028      10      FORMAT (' OLD, NEW CENTERS '2F10.2)
0029      TEST=CENA*0.25
0030      IF(ABS(DELNC).GT.TEST) GO TO 20
0031      IF(NIT.GT.10) GO TO 20
0032      NCENT=CEN+0.5
0033      K=NCENT
0034      PHD=CEN-FLOAT(NCENT)
0035      WRITE(5,169) AM0,BM0,DELNC,CENA,CEN,LIM
0036      169      FORMAT(2F20.5,3F7.3,I10)
0037      WRITE(6,169) AM0,BM0,DELNC,CENA,CEN,LIM
0038      IF((ABS(DELNC).LT.0.05).AND.(NIT.GT.1)) RETURN
0039      GO TO 1
0040      20      PRINT 11
0041      11      FORMAT (' SLOW CONVERGENCE IN COMBINE!')
0042      CEN=CENA+0.5*DELNC
0043      K=CEN+0.5
0044      PHD=CEN-FLOAT(K)
0045      RETURN
0046      END

```

ROUTINES CALLED:

IABS , MIN0 , FLOAT , COSTR , ABS

SWITCHES = /ON

58

BLOCK LENGTH

```

0001      FUNCTION DSUM (N,M,DATA)
)      C    T. TARBELL    JULY 12, 1971.
0002          DIMENSION DATA(512)
0003          SUM=0.0
0004          IF(M.LT.N) GO TO 2
0005          IF ((N.LE.0).OR.(N.GT.512)) GO TO 2
0006          IF ((M.LE.0).OR.(M.GT.512)) GO TO 2
0007          DO 1 I=N,M
0008      1      SUM=SUM+DATA(I)
0009              DSUM=SUM
0010          RETURN
0011      2      DSUM=0.
0012          RETURN
0013          END

```

SWITCHES = /ON

BLOCK	LENGTH
DSUM	140 (000430)*

```

**COMPILER --- CORE**
      PHASE      USED  FREE
DECLARATIVES  00366 13425
EXECUTABLES   00527 13264
ASSEMBLY      00988 15720

```

```

0001      SUBROUTINE ZER(DAT,N,NC,LIMIT,SLOP,WNEW,PHD,IZERO,NCOS)
C      ZER RETURNS THE FIRST POSITIVE ZERO OF THE SINE TRANSFORM OF DAT
C      USES METHOD OF FALSE POSITION WITH 5 ITERATIONS
C      WNEW IS VALUE OF TRANS. VARIABLE AT ZERO. IZERO=0 IF NO ZERO FOUND
C      ORIGIN IS AT NC, ENDPOINTS NC+LIMIT,NC-LIMIT. LENGTH IS N
C      SLOPE IS THE SIN TRANS. AT W = .1/LIMIT PHD IS PHASE FACTOR
C      IF NCOS LE 0, FINDS 1ST 2 ROOTS OF COS TRANS,
C      SLOPE=SMALLER ROOT, WNEW = LARGER ROOT
C      T. TARBELL JULY 22, 1971
0002      DIMENSION DAT(N)
0003      SLOP=0.
0004      WNEW=0.
0005      IF(NCOS.LE.0) GO TO 41
0006      W2=0.1/FLOAT(LIMIT)
0007      A=COSTR(DAT,N,NC,LIMIT,W2,SLOP,PHD,NCOS)
0008      SLOP=SLOP/(A*W2)
0009      41      CONTINUE
0010      W2=0.075
0011      TEM1=TRAN1(DAT,N,NC,LIMIT,W2,PHD,NCOS)
0012      W1=W2
0013      S1=SIGN(1.,TEM1)
0014      TEM2=0.
0015      1      W2=W2+0.025
0016      IF(W2.GE.1.) GO TO 4
0017      TEM2=TRAN1(DAT,N,NC,LIMIT,W2,PHD,NCOS)
0018      S2=SIGN(1.,TEM2)
0019      IF(S1.NE.S2) GO TO 20
0020      TEM1=TEM2
0021      S1=S2
0022      W1=W2
0023      GO TO 1
0024      4      CONTINUE
0025      PRINT 100
0026      100    FORMAT (' NO ZERO FOUND')
0027      IZERO=0
0028      RETURN
0029      20      CONTINUE
0030      TNEW=0.
0031      ITER=0
0032      21      ITER=ITER+1
0033      WNEW=W1-(W2-W1)*TEM1/(TEM2-TEM1)
0034      IF(ITER.GT.4) GO TO 30
0035      TNEW=TRAN1(DAT,N,NC,LIMIT,WNEW,PHD,NCOS)
0036      SNEW=SIGN(1.,TNEW)
0037      IF(SNEW.EQ.S1) GO TO 22
0038      S2=SNEW
0039      TEM2=TNEW
0040      W2=WNEW
0041      GO TO 21
0042      22      S1=SNEW
0043      TEM1=TNEW
0044      W1=WNEW
0045      GO TO 21
0046      30      CONTINUE
0047      IF(NCOS.GT.0) RETURN

```

```
0048      IF(SLOP*NE*0.) RETURN
0049      SLOP=WNEW
0050      W1=AMAX1(W1,W2)
0051      IF(W1*NE*W2) GO TO 43
0052      S1=S2
0053      TEM1=TEM2
0054      GO TO 1
0055 43      W2=W1
0056      GO TO 1
0057      END
```

ROUTINES CALLED:

FLOAT , COSTR , TRAN1 , SIGN , AMAX1

SWITCHES = /ON

BLOCK	LENGTH
ZER	615 (002316)*

```
**COMPILER --- CORE**
  PHASE      USED  FREE
DECLARATIVES 00366 13425
EXECUTABLES  00687 13104
ASSEMBLY     01364 15344
```

```

0001      FUNCTION COSTR(DATA,N,NC,LIMIT,W,SINTR,PHD,NCOS)
C      RETURNS THE COSINE TRANSFORM (COSTR) AND SINE TRANSFORM (SINTR)
C      USES SIN COS RECURSION FORMULAS FOR SPEED. NOT AN FFT. DONT INVERT
C      DATA IS THE ARRAY OF DIMENSION N, WITH ORIGIN AT NC AND ENDPOINTS
C      AT +/-LIMIT. W IS VALUE OF TRANSFORM VARIABLE. +SIN(W*X) IS USED
C      THE TRANSFORM IS THEN SHIFTED IN PHASE BY *PHD
C      T. TARBELL JULY 22, 1971
C      IF NCOS NEG., COSTR=SIN. TRAN. SINTR=COS. TRAN.
0002      DIMENSION DATA(N)
0003      CT=DATA(NC)
0004      IF(W.EQ.0.) GO TO 10
0005      ST=0.
0006      CTEM=1.0
0007      STEM=0.
0008      COSW=COS(W)
0009      COSW2=2.*COSW
0010      SINW=SIN(W)
0011      DO 1 I=1,LIMIT
0012      J=NC+I
0013      K=NC-I
0014      CT=CT+COSW*(DATA(J)+DATA(K))
0015      ST=ST+SINW*(DATA(J)-DATA(K))
0016      H=COSW2*COSW-CTEM
0017      CTEM=COSW
0018      COSW=H
0019      H=COSW2*SINW-STEM
0020      STEM=SINW
0021      3 SINW=H
0022      1 CONTINUE
0023      IF(PHD.NE.0.) GO TO 2
0024      SINTR=ST
0025      COSTR=CT
0026      IF(NCOS.GE.0) RETURN
0027      T=SINTR
0028      SINTR=COSTR
0029      COSTR=T
0030      RETURN
0031      2 PH=PHD*W
0032      CPH=COS(PH)
0033      SPH=SIN(PH)
0034      SINTR=CPH*ST+SPH*CT
0035      COSTR=CPH*CT+SPH*ST
0036      RETURN
0037      10 DO 11 I=1,LIMIT
0038      11 CT=(CT+DATA(NC+I)+DATA(NC-I))
0039      COSTR=CT
0040      SINTR=0.
0041      IF(NCOS.GE.0) RETURN
0042      SINTR=COSTR
0043      COSTR=0.
0044      RETURN
0045      END

```

ROUTINES CALLED:
COS , SIN

SWITCHES = /ON

BLOCK LENGTH


```

0001      FUNCTION TRAN1(DATA,N,NC,LIMIT,W,PHD,NCOS)
C      SAME AS COSTR, BUT COMPUTES ONLY ONE TRANSFORM
C      TRAN1 = SIN. TRAN. IF NCOS.GT.0
C      = COS. TRAN. IF NCOS.LE.0
C      TDT 4/5/74
0002      DIMENSION DATA(N)
0003      IF(NCOS.LE.0) GO TO 200
0004      PH=W*PHD
0005      STEM=-SIN(PH)
0006      ST=DATA(NC)*STEM
0007      COSW=COS(W)
0008      SINW=SIN(W)
0009      COSW2=2.*COSW
0010      DO 100 I=1,LIMIT
0011      J=NC+I
0012      K=NC-I
0013      ST=ST+SINW*(DATA(J)-DATA(K))
0014      H=COSW2*SINW-STEM
0015      STEM=SINW
0016      SINW=H
0017      100 CONTINUE
0018      TRAN1=ST
0019      RETURN
0020      200 CONTINUE
C      COS. TRAN. PART NOT NEEDED YET
0021      TRAN1=0.
0022      RETURN
0023      END

```

ROUTINES CALLED:

SIN , COS

SWITCHES = /ON

BLOCK	LENGTH
TRAN1	263 (001016)*

```

**COMPILER --*- CORE**
      PHASE      USED  FREE
DECLARATIVES  00446 13345
EXECUTABLES   00607 13184
ASSEMBLY      01084 15624

```

```

0001      SUBROUTINE POLYAN(DATA,N,NC,LIMIT,PHD,KMODE)
C      OVERSEES FOURIER TRANSFORM AND LEG. POL.
C      ANALYSIS OF CIRCULAR POARIZATIONS
C      IF ABS(COSGAM) GT 0.7, ALSO USES SEARS RELS. ON B/(A-C)
0002      DIMENSION DATA(N),TX(9),GPT(5),PLCOF(5),GWT(5,5),SPT(5)
0003      COMMON/RESULTS/X,COSGAM,GAM,SLOPE,IZERO,NORAT,
1      PLCOF,A,B,C,ALOR
0004      DATA SPT/1.0,0.91390,0.67095,0.31490,0./
0005      DATA NGPT/4/,GPT/1.5708,1.1529,0.7355,0.3204,
1      0./, GWT/ .20897950, 0., -.52244898, 0., .70530612,
2      .19091503, .23244581, -.24144505, -.59022408, -.21301324,
3      .13985270, .31111541, .22712270, -.09098065, -.45841186,
4      .06474248, .18434281, .27554684, .32346645, .31877204,
5      5+0./
0006      DATA PI/3.14159265/,RAD/57.2958/
0007      ACOS(Y)=ASIN(SQRT(1.-Y*Y))
0008      IZERO=5
0009      CALL ZER(DATA,N,NC,LIMIT,SLOPE,WNEW,PHD,IZERO,1)
0010      IF(KMODE.EQ.1) SLOPE=-SLOPE
0011      IF(IZERO.EQ.0) RETURN
0012      X=PI/WNEW
0013      DO 10 M=1,NGPT
0014      TX(M)=RAT(DATA,N,NC,LIMIT,PHD,GPT(M),SPT(M),X,NORAT)
0015      IF(NORAT.NE.0) RETURN
0016      10 CONTINUE
0017      DO 11 M=2,NGPT
0018      CKX=PI-GPT(M)
0019      TX(M+NGPT+1)=RAT(DATA,N,NC,LIMIT,PHD,CKX,SPT(M)
1      ,X,NORAT)
0020      IF(NORAT.NE.0) RETURN
0021      11 CONTINUE
0022      CALL PLINT(TX,PLCOF,NGPT,GWT)
0023      FAC=1./(PLCOF(1)+PLCOF(2))
0024      B=PLCOF(1)*FAC
0025      A=0.5*(1.+PLCOF(2))*FAC
0026      C=A-FAC
0027      IF(ABS(PLCOF(2)).GT.1.) GO TO 110
0028      COSGAM=1.*SIGN(1.,PLCOF(2))
0029      GO TO 199
0030      110 COSGAM=PLCOF(2)+SIGN(1.,PLCOF(2))*SQRT(PLCOF(2)*PLCOF(2)
1      -1.)
0031      199 IF(KMODE.EQ.1) COSGAM=-COSGAM
0032      200 CONTINUE
0033      RETURN
0034      END

```

ROUTINES CALLED:

ASIN , SQRT , ZER , RAT , PLINT , ABS , SIGN

SWITCHES = /ON

BLOCK	LENGTH
POLYAN 623	(002336)*
RESULTS 28	(000070)

COMPILER --- CORE
 PHASE USED FREE
 DECLARATIVES 00366 13425

```

0001      FUNCTION RAT(DATA,N,NC,LIMIT,PHD,PHAS,SKX,X,NORAT)
C          RAT = SIN(KX) * COSINE TRANS(K) / SIN TRANS(K)
C          SHOULD LOOK LIKE C1 + C2 * COS(KX)
C          NORAT.NE.0 IF INFINITY ENCOUNTERED
0002      DIMENSION DATA(N)
0003      DATA SMALL/1.E+4/
0004      NORAT=0
0005      W=PHAS/X
0006      CT=COSTR(DATA,N,NC,LIMIT,W,ST,PHD,1)
0007      SMAL=SMALL*ABS(CT)
0008      IF(ABS(ST).GT.SMAL) GO TO 10
0009      NORAT=10
0010      RAT=1.
0011      RETURN
0012  10      CONTINUE
0013      RAT=CT*SKX/ST
0014      RETURN
0015      END

```

ROUTINES CALLED:
COSTR , ABS

SWITCHES = /ON

BLOCK	LENGTH
RAT	183 (000556)*

```

**COMPILER --- CORE**
      PHASE      USED  FREE
DECLARATIVES  00366 13425
EXECUTABLES   00607 13184
ASSEMBLY      01020 15688

```

```

0001      SUBROUTINE PLINT(TX,PLCOF,NGPT,GWT)
C      PLCOF(NLEG) = COEFF. OF (NLEG+1)TH LEGENDRE POLY. IN
C      EXPANSION OF TX(COSKX)
C      GAUSSIAN QUAD. 2*NGPT-1 PTS.
0002      DIMENSION TX(9),PLCOF(5),GWT(5,5)
0003      DO 10 NLEG=1,5
0004      SUM=0.
0005      DO 1 I=1,NGPT
0006      1      SUM=SUM+TX(I)*GWT(NLEG,I)
0007      SGN=1.
0008      IF((NLEG.EQ.2).OR.(NLEG.EQ.4)) SGN=-1.
0009      NL=NGPT+1
0010      NU=NGPT+NGPT-1
0011      DO 2 I=NL,NU
0012      2      SUM=SUM+SGN*TX(I)*GWT(NLEG,I-NGPT+1)
0013      10      PLCOF(NLEG)=SUM
0014      RETURN
0015      END

```

SWITCHES = /ON

BLOCK	LENGTH
PLINT	214 (000654)*

```

**COMPILER --- CORE**
      PHASE      USED  FREE
DECLARATIVES 00446 13345
EXECUTABLES  00527 13264
ASSEMBLY     01044 15664

```

```

0001      SUBROUTINE BGOOD(GOLD,B,GAMF)
0002      DIMENSION GOLD(2,8)
0003      DATA DXMAX,DMAX,RESMAX,AVMIN/.25,.15,.15,150./
0004      PERCEN(X,Y)=2.*ABS(X-Y)/(X+Y+1.E-4)
0005      ACOS(Y)=ASIN(SQRT(1.-Y*Y))
0006      DATA PI,RAD/3.1415927,57.2958/
0007      DM=PERCEN(GOLD(1,4),GOLD(2,4))
0008      IF(DM.GT.DMAX) GO TO 100
0009      AVMOM=.5*(GOLD(1,4)+GOLD(2,4))
0010      IF(ABS(AVMOM).LT.AVMIN) GO TO 101
0011      ILOOK=1
0012      IF(GOLD(2,7).GT.RESMAX) GO TO 50
0013      ILOOK=2
0014      IF(GOLD(1,7).GT.RESMAX) GO TO 60
0015      GO TO 51
0016      50 IF(GOLD(1,7).LT.RESMAX) GO TO 60
0017      51 W1=1./PERCEN(GOLD(1,4),GOLD(1,5))
0018      W2=1./PERCEN(GOLD(2,4),GOLD(2,5))
0019      T=1./(W1+W2)
0020      W1=W1*T
0021      W2=W2*T
0022      B=W1*GOLD(1,2)+W2*GOLD(2,2)
0023      CGAMF=W1*GOLD(1,3)+W2*GOLD(2,3)
0024      GAMF=ACOS(CGAMF)
0025      IF(CGAMF.LE.0.) GAMF=PI-GAMF
0026      GAMF=GAMF*RAD
0027      RETURN
0028      60 B=GOLD(ILOOK,2)
0029      GAMF=ACOS(GOLD(ILOOK,3))
0030      IF(GOLD(ILOOK,3).LE.0.) GAMF=PI-GAMF
0031      GAMF=GAMF*RAD
0032      RETURN
0033      100 B=0.
0034      GAMF=0.
0035      RETURN
0036      101 B=ABS(AVMOM)
0037      GAMF=0.
0038      IF(AVMOM.LE.0.) GAMF=180.
0039      RETURN
0040      END

```

ROUTINES CALLED:

ABS , ASIN , SQRT

SWITCHES = /ON

BLOCK	LENGTH
BGOOD	594 (002244)*

```

**COMPILER --*-+ CORE**
  PHASE      USED  FREE
DECLARATIVES 00366 13425
EXECUTABLES  00687 13104
ASSEMBLY     01264 15444

```

```

C      PROGRAM MGRAM
C      READS PASS3 OUTPUT TAPE AND PRINTS MAPS OF RESULTS
C      NEEDS SUBROUTINES PACK, PMAPS
C      HANDLES ANY NO. PTS. ON SUN, PRINTED OUT IN BLOCKS
C      OF NMAX...TO CHANGE NMAX, CHANGE DIM.,EQ. & DATA
C      STATEMENTS FOR R, B, V, GAM, AND NMAX
C      TDT APRIL 5, 1974
0001      DIMENSION IFIL(5),R(4500),B(1500),GAM(1500),V(1500),
1          TITLE(3,5)
0002      COMMON IDATA(100),MAPS(5),DELX,DELY
0003      EQUIVALENCE (B(1),R(1)),(GAM(1),R(1501)),(V(1),
1          R(3001)),(XSTEP,IDATA(1)),(YSTEP,IDATA(3)),(NX,
2          IDATA(5)),(NY,IDATA(6))
0004      DATA R/4500*0./,NMAX/1500/
0005      DATA IBLANK/, /
0006      DATA TITLE/'B ',2*,' ', 'GA', 'M ', ' ', 'V ',2*,' ',
1          'CO', 'NT', ' ', 'BV', 'ER', 'T ' /
0007      5 WRITE (6,1010)
0008      1010 FORMAT(' TYPE FILENAME .CR TO ABORT!/')
0009      READ (6,1011) IFIL
0010      1011 FORMAT (5A2)
0011      IF(IFIL(1).EQ.IBLANK) GO TO 200
0012      WRITE (6,700)
0013      700 FORMAT(/' WHICH MAPS SHALL MOUSE PRINT?'/
1          ' TYPE 1 FOR YES, 0 FOR NO (I1 FORMAT)!/')
0014      DO 710 I=1,5
0015      WRITE (6,701) (TITLE(J,I),J=1,3)
0016      701 FORMAT(/' ',3A2/)
0017      READ(6,702) MAPS(I)
0018      702 FORMAT(I1)
0019      710 CONTINUE
0020      CALL PFEED(0)
0021      PRINT 1012,IFIL
0022      1012 FORMAT ('/' FILENAME ',.5A2//')
0023      CALL SETFIL(8,IFIL)
0024      1 READ(8,END=1000) ICODE,M3,(IDATA(M13),M13=1,M3)
0025      2 ICODE=ICODE+1
0026      GO TO (50,10,20) ICODE
0027      10 PRINT 1000,(IDATA(M),M=1,M3)
0028      WRITE(6,1000) (IDATA(M),M=1,M3)
0029      1000 FORMAT (500A2)
0030      GO TO 1
0031      20 NXTOT=NX
0032      NYTOT=NY
0033      DELX=XSTEP
0034      DELY=YSTEP
0035      21 IOVFLO=0
0036      IF(NXTOT*NYTOT.LE.NMAX) GO TO 22
0037      NY1ST=NYTOT
0038      IOVFLO=1
0039      NYTOT=INT(NMAX/NXTOT)
0040      22 CONTINUE
0041      PRINT 1102,NXTOT,NYTOT,DELX,DELY
0042      1102 FORMAT(/' NX, NY, DELX, DELY'2I5,2F7.3/)
0043      GO TO 1

```

```

0044      50      CALL PACK(R,NXTOT,NYTOT,B,GAM,V)
0045      IF(10VFLO.EQ.0) GO TO 5
0046      NYTOT=NY1ST=NYTOT
0047      GO TO 21
0048      100      END FILE 8
0049      PRINT 1004
0050      WRITE(6,1004)
0051      1004      FORMAT(/' EOF FOUND BY MAIN PROGRAM.'/)
0052      GO TO 5
0053      200      PRINT 1013
0054      WRITE(6,1013)
0055      1013      FORMAT(/' RUN ENDED!')
0056      STOP
0057      END

```

ROUTINES CALLED:

PFEED , SETFIL, INT , PACK

SWITCHES = /ON

BLOCK	LENGTH
MAIN. 9738	(046024)*
.\$\$\$\$ 109	(000332)

```

**COMPILER --- CORE**
      PHASE      USED  FREE
DECLARATIVES 00683 13108
EXECUTABLES  00767 13024
ASSEMBLY     01593 15115

```

```

0001      SUBROUTINE PACK(R,NX,NY,B,GAM,V)
0002      COMMON IDATA(100),MAPS(5),DELX,DELY
0003      DIMENSION R(2),JDATA(10),B(NX,NY),V(NX,NY),
1      GAM(NX,NY)
0004      EQUIVALENCE (BEE,IDATA(5)),(GAMA,IDATA(7)),
1      (VEE,IDATA(3)),(IX,IDATA(1)),(IY,IDATA(2))
0005      EQUIVALENCE (BE,JDATA(5)),(GAMM,JDATA(1)),
1      (VE,JDATA(3)),(JX,JDATA(1)),(JY,JDATA(2))
0006      B(IX,IY)=BEE
0007      GAM(IX,IY)=GAMA
0008      V(IX,IY)=VEE
0009      1      READ (R,END=100) ICODE,M3,(JDATA(M33),M33=1,M3)
0010      IF(ICODE.NE.0) GO TO 200
0011      B(JX,JY)=BE
0012      V(JX,JY)=VE
0013      GAM(JX,JY)=GAMM
0014      IF((JX.NE.NX).OR.(JY.NE.NY)) GO TO 1
0015      CALL PMAPS(B,GAM,V,NX,NY)
0016      RETURN
0017      100    PRINT 1000
0018      WRITE (6,1000)
0019      1000    FORMAT(/' PACK FOUND EOF'/)
0020      101    B(JX,JY)=BE
0021      V(JX,JY)=VE
0022      GAM(JX,JY)=GAMM
0023      CALL PMAPS(B,GAM,V,NX,JY)
0024      RETURN
0025      200    PRINT 1001, ICODE
0026      WRITE (6,1001) ICODE
0027      1001    FORMAT(/' DATA EXPECTED.,ICODE ='I5'/)
0028      IF(ICODE.EQ.1) WRITE (6,1002) JDATA
0029      IF(ICODE.EQ.1) PRINT 1002,JDATA
0030      1002    FORMAT(/' 40A2'/)
0031      CALL PMAPS(B,GAM,V,NX,JY)
0032      RETURN
0033      END

```

ROUTINES CALLED:
PMAPS

SWITCHES = /ON

BLOCK	LENGTH
PACK	485 (001712)*
.\$\$. 109	(000332)

COMPILER --- CORE

PHASE	USED	FREE
DECLARATIVES	00366	13425
EXECUTABLES	00793	12998
ASSEMBLY	01361	15347


```

0001 SUBROUTINE PMAPS(B,GAM,V,NX,NY)
0002 COMMON IDATA(100),MAPS(5),DELX,DELY
0003 DIMENSION SPACE(30)
0004 EQUIVALENCE (SPACE(1),IDATA(1))
0005 DIMENSION B(NX,NY),V(NX,NY),GAM(NX,NY)
0006 DATA RAD/57.2958/
0007 DO 100 I=1,5
0008 IF(MAPS(I)) 1,100,1
0009 1 GO TO (10,20,30,40,50) I
0010 10 CALL PFEED(0)
0011 PRINT 101,DELX,DELY
0012 101 FORMAT(/' MAG. FIELD (GAUSS) ON 'F7.3' BY '
1 F7.3' ARCSECOND GRID'//)
0013 DO 11 K=1,NY
0014 11 PRINT 102,(B(L,K),L=1,NX)
0015 102 FORMAT(/25F5.0)
0016 GO TO 100
0017 20 CALL PFEED(0)
0018 PRINT 201,DELX,DELY
0019 201 FORMAT(/' GAMMA (DEGREES) ON 'F7.3' BY 'F7.3
1 ' ARCSECOND GRID'//)
0020 DO 21 K=1,NY
0021 21 PRINT 102,(GAM(L,K),L=1,NX)
0022 GO TO 100
0023 30 CALL PFEED(0)
0024 PRINT 301,DELX,DELY
0025 301 FORMAT(/' REL. VEL. (KM/SEC: + = BLUESHIFT) ON 'F7.3
1 ' BY 'F7.3' ARCSECOND GRID'//)
0026 DO 31 K=1,NY
0027 31 PRINT 302,(V(L,K),L=1,NX)
0028 302 FORMAT(/25F5.1)
0029 GO TO 100
0030 40 CALL PFEED(0)
0031 PRINT 401,DELX,DELY
0032 401 FORMAT(/' CONTINUUM INTENSITY ON 'F7.3' BY '
1 F7.3' ARCSECOND GRID'//)
0033 DO 41 K=1,NY
0034 DO 42 L=1,NX
0035 42 SPACE(L)=100.*(B(L,K)/10.-INT(B(L,K)/10.))
0036 PRINT 102,(SPACE(L),L=1,NX)
0037 41 CONTINUE
0038 GO TO 100
0039 50 CALL PFEED(0)
0040 PRINT 501,DELX,DELY
0041 501 FORMAT(/' VERT. MAG. FIELD (GAUSS) ON 'F7.3
1 ' BY 'F7.3' ARCSECOND GRID'//)
0042 DO 51 K=1,NY
0043 DO 52 L=1,NX
0044 52 SPACE(L)=B(L,K)*COS(GAM(L,K)/RAD)
0045 PRINT 102,(SPACE(L),L=1,NX)
0046 51 CONTINUE
0047 100 CONTINUE
0048 RETURN
0049 END

```

ROUTINES CALLED:
PFEED , INT , COS